COMPONENTISTICA PROTESICA

ISTRUZIONI D'USO

IT

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Codice

PROSTHETIC COMPONENTS INSTRUCTIONS FOR USE

EN

INSTRUCTIONS FOR USE PROSTHETIC COMPONENTS

Prosthetic components consist of supporting devices that are connected to the implant; these devices are designed and constructed to be placed on dental implants so as to create an overall implant-prosthetic system capable of replacing one or more missing teeth. Fabricated prostheses are intended to restore correct masticatory, aesthetic and phonetic function to the user patients.

Prosthetic components are supplied in non-sterile pouches, sterilisation before use in the oral cavity is mandatory.

These operating instructions are provided in electronic form, which can be viewed at www.leadermedica.com.

Description

1. PROTESTIC ABUTMENTS

These are preformed accessories designed to support a ceramic or resin prosthesis. These abutments are applied to the biphasic implant during surgical reopening after the osseointegration phase. The straight abutment has a straight axis in relation to the axis of the implant. All abutments are attached to the biphasic implant by means of a connection screw. The abutments and connection screws are supplied in non-sterile pouches. They must be sterilised before use in the oral cavity.

- Straight and angled abutments: titanium abutments to which the prosthesis is attached
- Millable abutments: titanium abutments to which the prosthesis can be milled
- Locator abutments: titanium abutments for overdentures
- Temporary abutments: titanium abutments that allow the application of a temporary crown
- MUA: The MUA is an abutment. It represents a modality to ensure total retention and stability of a prosthesis, while the support remains divided between implants and tissue support. It consists of a transmucosal path, which can have different heights, and a cone representing the male of the attachment.

2. T-BASE/ PREMILLED/ CANNULA

The T-Base and Premilled are a base for crowns. The bases can also be cannulated.

3. PROTESTIVE SCREWS

These are the screws required to screw in abutments and superstructures. They are sold either as part of the abutment or separately as spare parts. Since they are used with abutments and superstructures, they are classified with the same class as the reference device.

4. HEALING SCREWS

The healing screw is a component that replaces the cap screw and is retained until the next prosthetic phase, allowing the gingiva to heal while preventing tissue regrowth within the implant. The maximum period of retention in the body is 180 days.

BOTM healing screws are made of methacrylic resins by 3D printing.

PRIMARY MATERIALS

MEDICAL DEVICES	PRIMARY MATERIALS
Prosthetic abutments	TITANIUM GR.5
T-base	TITANIUM GR.5
Bases	TITANIUM GR.5 / COBALT CHROME
Premilled	TITANIUM GR.5
Cannula	PMMA
Prosthetic screws	TITANIUM GR.5
Healing screws	TITANIUM GR.5
BOTM Healing screws	METHACRYLIC RESINS

CHEMICAL COMPOSITION OF MATERIAL

- Titanium Grade 5 Ti6Al4V - ASTM F136 and ISO 5832-3:

ELEMENT	MAXIMUM COMPOSITION LIMITS MASS FRACTION %
Aluminium	5.5 TO 6.75
Vanadium	3.5 TO 4.5
Iron	MAXIMUM 0.3
Oxygen	MAXIMUM 0.2
Carbon	MAXIMUM 0.08
Nitrogen	MAXIMUM 0.05
Hydrogen	MAXIMUM 0.015
Titanium	EQUILIBRIUM

- Methacrylic resin SprintRay EU Temporary Crown & Teeth:

ELEMENT	MAXIMUM COMPOSITION LIMITS MASS FRACTION %
Esterification products of 4,4'-isopropylidenediphenol, ethoxylated and 2-methylprop-2-enoic acid	25-45 %
7,7,9 (or 7,9,9)-trimethyl-4,13-dioxo-3,14-dioxo-5,12-diazahexadecane-1,16-diyl bismethacrylate	25-45 %
3,6,9-trioxaundecamethylenedimethacrylate	07-10 %
2-propenoic acid, reaction products with pentaerythritol	< 0.2 %
Diphenyl (2,4,6-trimethylbenzoyl)phosphine oxide	< 0.2 %

- COBALT CHROME:

ELEMENT	MAXIMUM COMPOSITION LIMITS MASS FRACTION %
Iron	MAXIMUM 0.75
Carbon	MAXIMUM 0.14
Nitrogen	MAXIMUM 0.25
Silicon	MAXIMUM 1.00
Manganese	MAXIMUM 1.00
Chromium	26.00 TO 30.00
Molybdenum	5.00 TO 7.00
Nickel	MAXIMUM 0.1
Cobalt	EQUILIBRIUM

Durability of the device

The average lifetime of prosthetic components is generally 10 years, unless undersized medical devices are used in relation to the implant site. In this case one has the obligation to inform the patient that for various reasons one had to operate outside the indications in the surgical protocol and thus of the possible risks to which the implant is subjected. This is one of the conditions under which the implant may have a shorter duration than normal. In addition, the patient must be thoroughly informed of the procedures to be followed for proper implant preservation. It is the doctor's duty to schedule routine check-ups and, if necessary, to carry out repairs on the implants to ensure that they function properly.

General risks

The use of devices is reserved for professionals who have acquired adequate professional training. Before using a component, the professional must check the mechanical integrity of the device and its packaging. In case of doubt about the use of the product, the professional should consult the local representative or the manufacturer. The manufacturer cannot be held responsible for modifications made to the device without its approval.

Absolute general contraindications

- Severe congenital or acquired mental or nervous illnesses, such that treatment and prognosis are problematic or impracticable.
- Recurrent or chronic nervous diseases in which the subject's relationship with reality goes through periods of lability, during which there is a possibility that the patient's psychic relationship with the endosseous prosthesis will be compromised.
- Serious cancer diseases with an inauspicious outcome.
- Serious diseases of: Bones-Connective system-Heart-Blood and Haemopoietic system-Endocrine system-Kidney-Liver-Nervous System.

Relative general contraindications

- Pregnancy and post-partum period: it is necessary to wait until the end of the pregnancy and lactation period.
- Analyses in the case of osteoporosis should be further investigated (acid and alkaline phosphatase blood values).
- In the case of non-serious heart disease, one must ascertain the seriousness of the pathology by consulting the cardiologist, then ascertain whether anticoagulants are being taken and, if necessary, in the case of multiple applications of implants, one must follow the cardiologist's judgement. Appropriate antibiotic coverage must be arranged. Alternatively, only one implant per operation can be applied in order to reduce blood loss, using a suture around the neck of the implant for haemostatic purposes.
- Acute joint rheumatism (RAA): operations are performed following the (antibiotic) prophylaxis prescribed by the attending physician.
- Allergic diathesis: an allergy to Titanium Gr. 5 must be ascertained.
- Essential trigeminal neuralgia: ascertain the severity of the syndrome, bearing in mind that the implant may prove to be a 'trigger', if the decision is made for the operation after prior authorisation by the treating neurologist.
- Anxiety-based motility alterations, clenching and grinding: stick to the concept of stress limitation. The use of ceramic as a material for prosthetic reconstruction is absolutely contraindicated.

Absolute local contraindications

Untreatable bone deficiency.

Relative local contraindications

- Anatomical defect correctable by reconstructive plastic.
- Acute or chronic local inflammation of soft tissue (parodontopathies).
- Acute or chronic local bone inflammatory processes.
- Presence of root debris at the surgery site.
- Presence of impacted teeth at the surgery site.

Warnings

1. Prosthetic components must only be inserted using 3P IMPLAFAVOURITE surgical instruments or Original lines.
 2. The original packaging of the 3P IMPLAFAVOURITE prosthetic components is a non-sterile bag with an identification label according to MDR.
 3. Prosthetic components are disposable and non-sterile and their reuse can lead to serious infections for the patient.
 4. Prosthetic components are intended to be sterilised before use. Wash in ultrasonic bath with detergent and disinfectant solution. Autoclave by steam cycle at 121°C for 15 minutes.
 5. The cycle used for sterilisation must be validated according to EN ISO 17665-1.
 6. The devices are intended for use in the patient placed immediately after cooling following sterilisation. The patient must be informed of the importance of cleaning and instructed in the hygiene and maintenance of the prosthesis. Patients with titanium dentures are prohibited from using toothpastes or mouthwashes containing free fluoride.
 7. If therapeutic irradiation of the head and neck is necessary, metal prostheses must be removed from the mouth.
 8. Devices as part of the prosthetic components were not evaluated for safety and compatibility in the MR environment. They were not tested for heating, migration or imaging artefacts in the MR environment. The safety of devices that are part of the prosthetic components in the MR environment is unknown. Scanning a patient with this device may cause injury to the patient.
 9. It is emphasised that CrCo alloy is a CMR material.
- In the event of sudden pain or complications, advise the patient to consult the doctor, surgeon or dentist without delay.

Storage
Store in a way that prevents damage to the packaging. Protect against moisture. Store at room temperature in a dry place. Do not reuse. Do not use the device if the packaging is damaged.

Disposal
The devices are used in a hospital environment and must therefore be disposed of according to the applicable waste disposal laws enforced by the hospital. In particular, used devices that may be polluted and biologically contaminated must be disposed of as special waste. Since they are intended for use in hospital facilities that normally implement such special collection, it is not necessary to provide such indications on labelling.

You must report any serious incidents occurring in connection with the medical device you supply to the manufacturer and to the competent authority of the Member State where you are located.

The safety and clinical performance summary is available at EUDAMED (<https://ec.europa.eu/Tools/eudamed>), identified by the basic UDI-DI codes, and until this platform is fully operational can be requested by sending an e-mail to info@leadermedica.com.

Symbols

	Caution!		Medical Device		Latex free
	Consult instructions for use		Date of manufacture		Keep dry
	Do not re-use		Manufactures		Device supplied non-sterile
	Batch code		Unique Device Identifier		
	Do not use if package is damaged		Code (Catalogue number)		

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