STRUMENTARIO CHIRURGICO

ISTRUZIONI D'USO

IT

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Codice

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Fabbricante

SURGICAL INSTRUMENTS INSTRUCTIONS FOR USE

EN

INSTRUCTIONS FOR USE SURGICAL INSTRUMENTS

Wrenches/Motor Wrenches, Micromotor Tap, Drills, Stoppers, Crestal Levelers, Extractors for Broken Screws, Mucotomes, Drill Extension: Medical device Reg. 2017/745 et seq., Class IIa, Annex IX.

3P IMPLAFAVOURITE instruments must only be used for 3P IMPLAFAVOURITE dental implants.

Lifetime

After 30 cycles of use, and the corresponding sterilisation cycles, the instruments may wear out to such an extent that they need to be replaced.

The service life of the instrument depends on the actual conditions of use and, in the case of unsuitable use, may be shorter.

In any case, after 100 uses, suitable qualities are no longer guaranteed and replacement is required.

General risks

The product is intended for use by appropriately trained medical professionals.

Before using the instrument, the user must check the condition of the cutting edges.

Please refer to surgical protocols for use. If in doubt about the use of the product, consult your local representative or the manufacturer.

The manufacturer cannot be held responsible for modifications made to the instrument without its approval.

Infection risks

Disinfect and sterilise before first use and after each subsequent use.

Use the reprocessing procedure described under "Sterilisation".

Electrical risks

The instruments may only be used in connection with handpieces connected to low-voltage units for medical use in accordance with the Medical Device Regulation MDR 2017/745.

Mechanical risks

Check the tightness of the instrument after each connection to the handpiece.

Do not insert or remove the instrument from a rotating motor. Never handle the rotating instrument during operation. Never operate the instrument exchange system during operation and before the rotation of the drill or tap has come to a complete stop.

The implantable device should never be used in an osteotomy excavated with a larger diameter drill.

Sterilisation

INSTRUCTIONS FOR REPROCESSING DENTAL INSTRUMENTS			
	PROCESS	PROCESS STEPS	DESCRIPTION, INSTRUCTIONS, WARNINGS
1.	Initial treatment at the point of use	1.1 Soaking	Instrument soaking should take place after use within a maximum of 2 hours. Use products with a proven enzyme component (e.g. DGHM, FDA approval, or CE mark) that are suitable for instrument-compatible disinfection. Use the disinfectant according to the product manufacturer's instructions. Use a synthetic cloth or soft brush for preliminary removal of dirt residues. Rinse the instrument under running water. WARNING: The disinfectant used for soaking is only for personal protection and does not replace the subsequent cleaning and disinfection steps.
		1.2 Containment and transport	No special requirements. It is recommended to process the instruments as soon as reasonably practical after use and soaking.
2.	Preparation before cleaning	-	No special requirements
3.	Automated cleaning, disinfection and drying	Thermal	A mechanical procedure using a washer-disinfector (washing, disinfecting and drying apparatus) is recommended for cleaning instruments. We recommend the use of neutral or slightly alkaline cleaning agents, free of critical substances (depending on concentration). If alkaline cleaning agents (pH 9.5 - 11.5) are used, colour changes on metal surfaces are possible. However, this does not impair the functionality of the instrument. Avoid strongly alkaline cleaning agents (pH > 11.5). The suitability of the products for effective cleaning and disinfection by means of a mechanical procedure (93°C, 10 min) using alkaline detergents and the addition of a surfactant agent has been verified. The use of other (or non-equivalent) cleaning agents is not the responsibility of the manufacturer. Only use cleaning agents proven to be effective and comparable to those presented in this section. When choosing the cleaning system, make sure that: - in principle, it is suitable for cleaning instruments; - the chemical agents are compatible with the instruments. Observe the concentrations and action times specified by the manufacturer of the cleaning agent. Immerse the unloaded instrument, with the lid open, in the cleaning and disinfection device and follow the instructions of the device manufacturer.
4.	Manual cleaning and disinfection	-	The use of products with enzyme and disinfectant components is recommended. Use the product according to the manufacturer's instructions. Immerse the instrument with the lid open downwards in the solution in the ultrasonic bath and use a soft brush to remove any residue remaining from the procedure. After treatment, rinse the instrument under running water. Check the cleanliness and, if dirt and/or foreign material are still visible, repeat the procedure.
5.	Manual drying	-	Dry the instrument thoroughly with alcohol-based disposable disinfectant wipes or a sterile, lint-free cloth and, if necessary, using compressed air. To avoid damaging the products, the air pressure must be < 3 bar. It is advisable that the air used for drying is filtered.
6.	Maintenance inspection and testing	-	Check the functionality of the products in accordance with the respective instructions for use and exclude visible damage. For moving or rotating parts, the use of instrument oil is not recommended, as certain plastics tend to expand and oils adversely affect proper functioning
7.	Packaging	-	Prior to sterilisation, the instrument must be placed in sterilisable single-use medical grade packaging (single or double packaging) and/or in a suitable sterilisation container: # compliant with EN 868/ ISO 11607; # suitable for steam sterilisation (temperature resistance up to 134°C, sufficient steam permeability); # regularly maintained (sterilisation container) CAUTION: Do not sterilise the instrument inside the original transport packaging.
8.	Sterilizzazione	In steam autoclave (Vacuum cycle at 134°C for 4 min)	Use only the following procedures for sterilisation; do not use any other procedures. Steam sterilisation, check that: # the steam steriliser (vacuum type) complies with EN 13060 and EN 285. # procedure is approved according to ISO 17665. # max. sterilisation temperature 134°C (plus tolerance according to ISO 17665) # sterilisation time: 4 minutes at 134°C. The suitability of the products for effective sterilisation with the above sterilisation times and temperatures has been verified. Do not use hot air sterilisation.
9.	Conservation	-	After sterilisation, store the products in a cool dry place and keep them in the sterilisation packaging.
10.	Checking and testing operation	-	Before each application, carry out a visual check and function test as per the operating instructions.

DEVICE DESCRIPTION

Surgical instruments dedicated to 3P IMPLAFAVOURITE implant systems are medical devices intended for use in the oral cavity, for temporary use (continuous duration not exceeding 60 minutes) and reusable after washing and sterilisation.

The surgical instruments are manufactured according to the most stringent European and American standards and are intended for use with dental implants manufactured by 3P IMPLAFAVOURITE.

DRILLS

Surgical drills are made of surgical steel characterised by high resistance to corrosion and wear. They are intended for use with a contra-angle handpiece attachment and must be used with a suitable micromotor.

Lever movements increase the risk of fracture of the drills and must therefore be avoided. Sudden changes in speed must be avoided and pressure must never be applied to forcibly stop the rotation of the instrument, as this could develop excessive heat in the tissue being cut and premature wear of the instrument and micromotor.

It is always recommended to work intermittently with an oscillating movement in a vertical direction, in order to avoid overheating and wear of the working part and high heat increase in the alveoli affected by the cutting action of the drill.

The use of appropriate coolant is recommended.

Lack of adequate irrigation can lead to bone necrosis.

Drill wear is largely dependent on the type and density of the bone being drilled: harder bone results in greater instrument wear.

It is important and necessary to check burs for wear every 30 uses.

Drills should never be resharpened before use. Never use damaged, bent or worn instruments.

Insert the pilot drill into the handpiece, screw in the stop corresponding to the height of the chosen implant and mill. The recommended rotation speed is 600 rpm (**never exceed 1000 rpm**).

Proceed with the preparation of the implant tunnel by progressively enlarging the hole as indicated in the surgical protocol, carefully checking the depth reached each time using the appropriate gauge. To avoid dangerous overheating in the bone, the rotation speed **must never exceed 1000 rpm**.

Time in the patient's mouth: temporary duration.

TAPPERS

These are cutting devices to be used in the case of very dense bone to mask the implant site.

Time in the patient's mouth: temporary duration.

HANDLE/ MICROMOTOR KEYS

These are devices powered by a micromotor, used to screw in implants quickly and constantly. The device does not come into contact with the body.

Time in the patient's mouth: temporary duration.

STOPPER

These are devices that allow precise depth control, preventing over-drilling and reducing the risk of injury to adjacent structures.

Time in the patient's mouth: temporary duration.

CRESTAL LEVELERS

These are devices that allow the bone crest to be levelled to allow better use of subsequent instruments.

Dwell time in the patient's mouth: temporary duration.

EXTRACTORS FOR BROKEN SCREWS

These are devices used to remove damaged or broken screws from dental implants. These instruments are essential to ensure the safety and effectiveness of dental implant maintenance and repair procedures.

Time in the patient's mouth: temporary duration.

MUCOTOMI

These are instruments used to create a complete incision of the flap of the muscularis.

Time in patient's mouth: temporary duration.

EXTENSION FOR DRILLS

This is a tool to be used if the total length of the drills is too short. The device does not come into contact with the body.

Dwell time in the patient's mouth: temporary duration

Warnings

It is recommended that surgical drills be used for a maximum of 100 performances, beyond which their suitable quality is no longer guaranteed.

To ensure that site preparation is carried out without excessive trauma to the bone and does not delay or compromise healing, it is necessary:

- Always bear in mind that burs overprepare in length.
- Bear in mind that all indications of drill depth lead to implant placement at a predetermined level during the design phase; if placement at a lower level is necessary, this must therefore be taken into account during all stages of site preparation.
- To maintain a low bone temperature and avoid necrosis of the insulted cells, provide good external irrigation of the drills and tappers, using saline cooled to approximately 4°C.
- Use motors with controllable RPM and power. Do not exert constant pressure on the handpiece up to the maximum depth, but drill with a piston rotation, advancing and retracting a few millimetres at a time, in order to favour the inflow of the saline solution and allow the evacuation of residues deriving from the drilling and tapping of the bone.- If abnormal resistance is felt before reaching the planned depth, it is a good idea to stop, check direction and depth and, if necessary, take an X-ray with the drill inserted to compare it with the reference overview used for the implant design.
- Use intact and sharp drills and taps.
- Keep the flap well spread and avoid contact with the rotary drills.
- Ensure good suction of blood and excess saline to keep the area in which you are operating well cleansed and visible.
- Control the thrust on the handpiece well to prevent sudden failure from causing serious damage to nerve or vascular structures.
- Only to be used with 3P IMPLAFAVOURITE components and their original lines. the saline solution and allow the evacuation of the residues deriving from the drilling and tapping of t

STORAGE

Store in a way that prevents damage to the packaging. Protect against moisture. Store at room temperature in a dry place. Do not use the device if the packaging is damaged.

DISPOSAL

The devices are used in a hospital environment and must therefore be disposed of according to the applicable waste disposal laws enforced by the hospital. In particular, used devices that may be polluted and biologically contaminated must be disposed of as special waste. As they are intended to be used in hospital facilities that normally implement such special collection, it is not necessary to provide such indications on labelling.

Limitation of Liability: 3P IMPLAFAVOURITE Srl shall not be liable for any direct, indirect and/or incidental damages arising from the application of the device, whether resulting from failure to check functionality, sterilisation, validity, surgical intervention or failure to comply with the above Instructions for Use.

You must report any serious incident occurring in connection with the medical device supplied by us to the manufacturer and to the competent authority of the Member State in which you are located.

Symbols

	Caution!		Do not use if package is damaged and consult instructions for use		Unique Device Identifier
	Consult instructions for use		Medical Device		Code (Catalogue number)
	Batch code		Date of manufacture		Keep dry
	Non-Sterile		Manufactures		

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