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## *Surgical protocol, insertion sequence*

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*The principles of B.O.P.T.<sup>®</sup> have been successfully extended to implantology, creating a specific surgical sequence for the IL<sup>®</sup> and IL WIDE<sup>®</sup> implant lines. IL<sup>®</sup> implant placement guidelines take into account soft tissue thickness and the anatomical conformation of the bone crest. This integrated approach allows for a prosthetically guided relationship in implant placement. The formation of circular collagen fibers in the microgrooves of the convergent neck, along with the formation of a short epithelial attachment around the prosthetic device, guarantee the long-term stability of the IL<sup>®</sup> implant.*

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*For illustrative purposes, the surgical sequence for an implant with the following characteristics is shown below:*

**Ø 3.80 MM, H 10.8 mm, NECK 2.8 mm - STOP USED  
FOR THE 10 mm HIGH IMPLANT.**

	PRIMA FRESA FIRST DRILL	SECONDA FRESA SECOND DRILL	TERZA FRESA THIRD DRILL	QUARTA FRESA FOUR DRILL	QUINTA FRESA FIFTH DRILL
Ø EFFETTIVO FRESA ACTUAL DRILL	LANCEOLATA LANCEOLATE	2.0	2.5	2.8	3.0
Ø IMPIANTO IMPLANT	-	-	-	3.4	3.8

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### IMPLANT INSERTION

*The IL<sup>®</sup> implant is a tissue-level implant and therefore must be placed with the convergent neck portion outside the bone. This maximizes the advantages of having the abutment/prosthesis connection far from the bone and the collagen fibers can create a biological seal around the implant neck.*

*The clinician decides on the placement based on the quantity and quality of the soft tissue, taking care to leave the connection slightly below the gingival margin.*

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### NECK CHOICE

*The bone crest and gingival thickness can vary when choosing between intramucosal or tissue-level insertion that is why we have introduced four neck heights: 0.9, 1.8, 2.8, and 3.8, to adapt to any situation.*

*The implant can be placed in a tissue-level position with minimal controlled soft tissue invasion, or in a slightly deeper intramucosal position. The decision will be made by the clinician during the case planning phase, also taking into account the prosthesis with which the patient will be rehabilitated (i.e. cemented, screw-retained, single or multiple implants). Correct implant placement should take into account the CEJ of the adjacent tooth, if it is present, aligning the connection with this line. In the case of more extensive edentulism, the location of the new crown must be considered so as to respect the concept of biological width (the minimum mandatory requirement is approximately 2 mm).*

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### USE OF THE STOP AND THE COMPENSATION SPACE

*The use of the stop is designed to allow the first 0.8 mm of the neck to sink into the bone.*

*This surface, called the compensation space, was designed to adapt to and compensate for the inevitable post-surgical bone remodeling.*

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### THE IL WIDE<sup>®</sup> IMPLANT



*For greater stability and mechanical strength in the posterior quadrants, we recommend using the IL WIDE<sup>®</sup> implant with the same insertion procedures as the IL<sup>®</sup> implant.*

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
## THE COMPENSATION SPACE

When measuring bone availability on radiographic images, consider that the drill action includes 0.8 mm, which must be added to the size of the chosen implant. The 0.8 mm refers to the compensation space provided by the protocol.

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To accurately determine the drilling depth, it is recommended to consider the length of the drill bit, which varies for each diameter.

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
*In keeping with the close synergy between the B.O.P.T.<sup>®</sup> philosophy and the convergent implantology of the IL<sup>®</sup> line, Leader Medica has developed a series of dedicated B.O.P.T.<sup>®</sup> preparation kits. These kits were designed in direct collaboration with Professor Ignazio Loi and manufactured according to his specific clinical indications and operating protocols with the aim of offering perfectly integrated and biologically oriented surgical and prosthetic solutions.*

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*The vertical preparation is created, using dedicated diamond burs of various shapes and grit, depending on the step of the preparation they are going to be used for.*


### PREPARATION AND FINISHING DRILLS FOR TEMPORARY BRIDGES

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*The sequential use of drills, discs, rubber, and brush allows for specific use in the various steps until the final, finished and polished temporary restoration is created. Each drill has a specific function in each step.*

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*The indications and images are just general instructions that should be adapted on a case-by-case basis according to the patient's needs. These recommendations are not intended to replace specialized training or the expertise of the clinician, who—also based on personal experience—can make different decisions.*

### **ESSENTIAL PROCEDURES TO ENSURE THE SAFETY, EFFECTIVENESS, AND FUNCTIONALITY OF INSTRUMENTS**

*All surgical instruments must be thoroughly cleaned, disinfected, and sterilized before first use and before each subsequent reuse. It is essential to check the proper functioning of the instruments before each procedure and, if necessary, to immediately replace any instruments showing signs of wear or malfunction. Failure to follow these instructions may increase the risk of infection and intraoperative complications for the patient.*

### **PREPARING THE AREA FOR IMPLANT PLACEMENT**

*Every procedure must be performed gradually and carefully, avoiding sudden movements that could compromise safety. For traditional surgery check the correct positioning of the drill stops, ensuring they are fully tightened. Choose the length in accordance with the implant insertion plan and the given warnings. Periodically check surgical micromotors used with drills and taps to prevent malfunctions. Regularly check the drills for wear and replace them if they lose effectiveness in order to ensure precision and safety during implant site preparation.*

### **PROSTHESIS MAINTENANCE**

*To prevent potential mechanical and biological complications associated with prosthetic implants—which are widely documented in the literature—it is essential for the patient to maintain proper oral hygiene and undergo scheduled periodic checkups. These measures help prolong the functional life of the device. Periodic checks of the calibration of abutment and prosthetic screws are also recommended to ensure the stability and efficiency of the device. If a feeling of instability or malfunction occurs, the patient should be advised to seek prompt medical attention.*

### **DISPOSING OF INSTRUMENTS AND IMPLANTS**

*Surgical instruments, consisting of small metal parts, can be disposed of as metal waste only if properly cleaned. Otherwise, they must be treated as biological waste. Similarly, dental implants removed from the patient's mouth must be treated as biological waste. For each disposal procedure, it is recommended to refer to the local regulations in force regarding the management of medical waste.*