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### *From analog to digital...*

*Precision, clinical predictability and aesthetics are at the core of our digital workflow. Each step of the process is carefully designed to effectively meet the needs of patients who require rapid and minimally invasive solutions with short recovery times, as well as professionals who seek maximum operational safety and reliable results, both functional and aesthetic.*

“ The meeting point between surgeon, dental technician, and patient ”

*Thanks to virtual planning integrated with 3D intraoral scanning, it is possible to acquire extremely accurate data, faithful to the patient's anatomy, and develop a personalized prosthetic design. Prosthetically guided surgery enables more controlled, predictable, and time-optimized procedures, improving the clinical experience and patient comfort.*

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*Data acquisition via intraoral scanner STL format and CT or CBCT in DICOM format.*

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*Implant and prosthetic design using 3PGuide software.*

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*3D design and production of the surgical guide, anatomical model, and temporary prosthesis.*

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*Prosthetically Guided Implant Surgery.*

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*With the 3P Guide software from Leader Medica, we provide you with a powerful tool for a comprehensive diagnosis and extremely accurate surgical and prosthetic planning, allowing you to complete the entire workflow.*

**Flexible:** Compatible with all DICOM formats, it easily adapts to any radiological examination. It is also compatible with Mac, Windows and iOS.

**Integrable:** it interfaces with major CAD/CAM systems for fluid and consistent modeling.

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*Interconnected modules for a seamless workflow*

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*Integrated cloud-based file sharing platform*

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*Intuitive procedures, thanks to automated functions*

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*Open architecture for maximum flexibility*

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*Integrated artificial intelligence*

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*Why choose surgically guided implantology?*

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01

## SAFETY DURING SURGICAL PROCESSING

*The combination of intraoral scanning and 3D radiology, integrated into digital planning, allows for an accurate analysis of sensitive anatomical structures (such as nerves, vessels, and bone structure), ensuring safer, more controlled, and highly predictable procedures.*

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02

## PREDICTABLE FINAL RESULTS

*Digital tools allow for precise visualization of the final prosthetic results already in the pre-operative phase, improving the planning and effectiveness of the entire treatment.*

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03

## MAXIMUM OPERATIONAL PRECISION

*The customized surgical guide, created with a 3D printer, follows the patient's anatomy, ensuring the correct implant placement and reducing the margin for error and post-surgical complications.*

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04

## SHARED CLOUD SOFTWARE

*The integration of digital platforms with cloud access facilitates seamless and immediate communication between the clinician and the laboratory, allowing for real-time project sharing.*

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

05

## DEDICATED TECHNICAL SUPPORT

*A team of specialists provides ongoing assistance for the use of digital systems.*

# Protocollo Chirurgia Protesicamente Guidata

## *Prosthetically Guided Surgery Protocol*

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### MEDICAL HISTORY AND DENTAL EXAMINATION

The medical history is used to collect all the patient's data to allow for a correct diagnosis. The dental examination aims to assess the general health of the mouth, preventative care, and assess the operating spaces, as a limited oral opening can impede the correct execution of the technique or limit the passage of drills. An opening of at least 40 mm is recommended.

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### DATA ACQUISITION - PARTIAL EDENTULOUS PROTOCOL

#### FILE ACQUISITION FROM AN INTRAORAL SCANNER.

The digital process involves the use of an intraoral scanner to improve precision and patient comfort. Digital impressions can only be taken in cases of partial edentulism, with the goal of obtaining a suitable oral record of the soft tissues and tooth surfaces.

Scan files must be saved in STL/PLY/OBJ format.

To achieve optimal scan quality, follow the specific protocol for the intraoral scanner used. Scans must be free of holes and well-extended over the fornix and palate in the case of the upper arch.

#### FILE ACQUISITION FROM A 3D SCANNER

Starting with an analog workflow, the initial plaster model is taken using high-precision material, ensuring full extension, including the palate and retro-molar tuberosities.

The resulting plaster model is then scanned using a benchtop scanner, with careful attention to proper alignment with the indexes and, if applicable, with the prosthesis, reference markers, or diagnostic wax-ups. The entire dataset must be exported properly aligned in .STL format.

#### DIAGNOSTIC WAX-UP (IF APPLICABLE)

The diagnostic wax-up provides information about the crowns to be restored.

Having a diagnostic wax-up available for the case will simplify planning the aesthetic aspects of the implants and abutments, allowing for correct guided prosthetic planning. It can be created directly in digital format or derived from a scan of an analog wax-up.

The required digital format for the wax-up is STL.

#### CT SCAN ACQUISITION

The CT scan is performed by acquiring images of the arch, taking care to place special instruments between the two arches to avoid occlusion.

#### CT OR CBCT - DICOM FILE

The volumetric scan provides information about the bone, teeth, and other anatomical structures.

When scanning a single arch, the maxilla or mandible and part of the opposing bone must be completely visible.

This will make aesthetic assessment easier during planning review. Making the ostiomeatal complex visible will also be useful for evaluating sinus pathologies, if any.

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#### VOLUMETRIC SCAN SETTINGS:

Scan images must be saved in DICOM format. To achieve optimal image quality, check the device's specific DICOM export protocol. If an export protocol for volumetric scanning is not available, it is recommended to save DICOM images with a recommended axial slice thickness of 0.2 or 0.4 mm.

### DATA ACQUISITION - TOTAL EDENTULOUS PROTOCOL WITH DOUBLE SCAN

#### PREPARATION OF THE SCAN PROSTHESIS

The scan prosthesis is a pre-existing removable prosthesis or a radiolucent duplicate of the diagnostic dental setup with eight markers inserted or glued to the surface (e.g., quartz spheres, Fig. 1-2). The scan prosthesis is used to visualize the desired final prosthetic goal and to design the surgical guide.

1. Verify the pre-existing scan denture. A pre-existing denture can be reused if three conditions are met:
  - a) The dental setup is the final goal;
  - b) There are no metal parts, or the metal parts are removable for scanning;
  - c) The denture fits perfectly or can be modified to perfectly fit in the patient's mouth.

Fig. 2

If the above criteria are not met, a radiolucent duplicate of the diagnostic dental setup will be needed. Specifically, a relined denture is required.

2. Create eight small, evenly distributed surface cavities on the denture (on the vertical, buccolingual, and mesiodistal portions).
3. Use a drop of radiolucent resin and dip a marker into it. Place the resin with the marker in one of the cavities. If necessary, light-cure to reduce the curing time. Repeat this step for all markers.

**IMPORTANT:** The scanning denture must be radiolucent and free of metal parts.

#### PREPARING THE OCCLUSAL INDEX

The occlusal index stabilizes and separates the maxilla and mandible during the patient scan. Scans showing maxillary and/or mandibular movement cannot be used to design the guide.

How to prepare an occlusal index:

1. Apply radiolucent material to the patient's occlusal surface.
2. The patient should close their mouth, but not completely. The occlusal index should be 3-5 mm thick; subsequently remove any excess material with a sharp instrument.

Caution! Do not discard the occlusal index after scanning; the occlusal index can be reused to accurately position the guide.

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### CT OR CBCT - DICOM FILE

The volumetric scan provides information about the bone, teeth, and other anatomical structures. In cases of complete edentulism, two scans are required, as per the double scan protocol.

How to Perform a Double Scan:

#### Scan 1

Scan the patient with the occlusal index and the scanning prosthesis with the markers.

If you wish to use an external marker, the marker must be assembled with the scanning prosthesis, always securing it with radiolucent material.

The maxilla or mandible, including the occlusal index and part of the opposing bone, must be completely visible. This will facilitate the aesthetic assessment during planning. Making the ostiomeatal complex visible will also be useful for evaluating sinus pathologies, if any.

#### Scan 2

Scan the scanning prosthesis alone in the same right/left, upper/lower position used in the first scan.

The support material used for the scanning prosthesis must be more radiolucent than the prosthesis itself. Polyethylene and polyurethane foam are types of materials that can be used when performing a cone beam scan.

### MATCHING

After importing the DICOM of the CT scan into the planning software, and before performing the implant planning, it is necessary to register the STLs according to the type of protocol.

#### PARTIAL EDENTULOUS PROTOCOL


CT - Arch STL (in the case of intraoral scanning or benchtop scanning)

#### DOUBLE SCAN PROTOCOL

CT (patient + prosthesis + markers) - CT (prosthesis + markers)

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The neck section (H) must be inserted into the bone for implants with a 1.8, 2.8 or 3.8 mm neck (N). Implants with a 0.9 mm neck (N) are bone-level implants.

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The length of the drills and the offset of the cannulae allow for guided implant insertion as outlined in this protocol.

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The design of the drills allows the implant to be inserted into the bone an additional 2 mm compared to the protocol for specific clinical needs.

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The use of the IL\* implant with a 3.4 mm diameter must be carefully evaluated by the surgeon. It is indicated exclusively for placement in aesthetic areas, characterized by limited occlusal loading and angulations close to the implant axis.

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## 3D DESIGN AND PRODUCTION OF THE GUIDE

Based on the planning, the surgical guide is designed. The software allows you to define parameters such as offset, thickness, and undercut management. When designing the guide, it is important to consider the presence and size of the bushings, their spacing, and any support areas that may affect stability or create pre-contact with any existing teeth. In cases of total edentulism, the use of lateral pins and their dedicated drills is mandatory. The surgical guide, designed in the software, is exported in STL format and printed with a 3D printer.

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## STANDARD DRILL SEQUENCES AND IMPLANT SURGERY

### ANESTHESIA

It is important not to create anesthetic boluses that could prevent the proper seating of the guide.

### MUCOTOMY

Performed using a special Tissue Punch after seating the guide.

### OPERCULECTOMY

The surgical guide must be removed to facilitate the removal of the mucosal operculum.

### CRESTAL PREPARATION AND INITIAL DRILL

It creates the crestal plate for the implant head by removing all interference related to the alveolar ridge. The IL® guided kit includes both a dedicated crestal mill for each type of neck and the initial drill.

### DRILLING SEQUENCE

Starting with the first 2x6 mm drill, the subsequent longer drills are guided by both the surgical guide bushings and the 6 mm pilot drill.

The drilling sequence may vary, depending on implant types and the bone quality at the implant site.


On the left, the different recommended sequences:

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### The Crestal Mill has no stop.

It features dual-purpose notch at 8mm and 9 mm. It is possible to reach 10 mm based on the coronal surface of the mill.

### IMPLANT PLACEMENT

Special guided mounts screwed into the implant are used to position the implant while maintaining the insertion axis as given by the software design. These mounts also allow you to track the position of the connection system. Different mounts are available, depending on the neck length of the IL® or IL WIDE® implant.

### REMOVAL OF MOUNTS AND FIXING SYSTEMS: PROSTHETIC PHASE

Removal of mounts and, if applicable, insert or remove healing screws or first-stage screws.

### IMMEDIATE LOADING

The decision to perform immediate loading is at the discretion of the professional performing the procedure.

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
## Single body

- 100% Stainless Steel
- Secure Closure
- Machine Washable Design
- Autoclave sterilisation at 134°C
- Removable rubber pads for easy cleaning
- Removable tray
- Tilting surface
- Lightweight

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## GUIDE SLEEVES

*The guide sleeves are hexagonal-headed cylinders (indexed) that are incorporated into the surgical guide. Their purpose is to guide the instruments during surgery. By ensuring a fixed ratio of 9.00 mm between the instrument contact point on the bushing and the implant portion to be inserted into the bone, the planning created with the 3P GUIDE software is respected.*

*The guide sleeves are 4.00 mm high and come in two different diameters. The LARGE sleeve has an internal diameter of 5.50 mm, while the SMALL sleeve has an internal diameter of 4.15 mm. The SMALL sleeve is only suitable for implants with a diameter of 3.40 mm or 3.80 mm.*

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













## PIN SLEEVES

*Pin sleeves are metal cylinders incorporated into the surgical guide and are used to stabilize the guide, especially in cases where the guide cannot support sufficient residual elements. The guide must be stable to ensure safety during surgery. The absence of micromovements allows the desired accuracy to be achieved.*

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SEQUENZA PASSAGGI / STEP SEQUENCE	H 6.0	H 8.5	H 10.0	H 11.5	H 13.0	H 15.0
	1. tissue punch	1. tissue punch	1. tissue punch	1. tissue punch	1. tissue punch	1. tissue punch
2.	crestal mill 3,40	crestal mill 3,40	crestal mill 3,40	crestal mill 3,40	crestal mill 3,40	crestal mill 3,40
3.	initial drill	initial drill	initial drill	initial drill	initial drill	initial drill
SEQUENZA FRESE / DRILLS SEQUENCE (Ø x H)	2x6	2x6	2x6	2x6	2x6	2x6
	3,4x6	2x8,5	2x8,5	2x8,5	2x8,5	2x8,5
		3,4x8,5	2x10	2x10	2x10	2x10
			3,4x10	2x11,5	2x11,5	2x11,5
				3,4x10	2x13	2x13
				3,4x11,5	3,4x10	2x15
					3,4x11,5	3,4x10
					3,4x13	3,4x11,5
						3,4x13
						3,4x15

Ø 3.80 mm 

MOUNTER COLLO / NECK IL® 1,8 - 2,8 - 3,8 / IL WIDE® 0,9 - 1,8 - 2,8  
TIPO DI BOCCOLA / SLEEVE TYPE 4,15 - 5,50

SEQUENZA PASSAGGI / STEP SEQUENCE	H 6.0	H 8.5	H 10.0	H 11.5	H 13.0	H 15.0
	1. tissue punch	1. tissue punch	1. tissue punch	1. tissue punch	1. tissue punch	1. tissue punch
2.	crestal mill 3,80	crestal mill 3,80	crestal mill 3,80	crestal mill 3,80	crestal mill 3,80	crestal mill 3,80
3.	initial drill	initial drill	initial drill	initial drill	initial drill	initial drill
SEQUENZA FRESE / DRILLS SEQUENCE (Ø x H)	2x6	2x6	2x6	2x6	2x6	2x6
	3,4x6	2x8,5	2x8,5	2x8,5	2x8,5	2x8,5
	3,8x6	3,4x8,5	2x10	2x10	2x10	2x10
		3,8x8,5	3,4x10	2x11,5	2x11,5	2x11,5
			3,8x10	3,4x10	2x13	2x13
				3,4x11,5	3,4x10	2x15
				3,8x10	3,4x11,5	3,4x10
				3,8x11,5	3,4x13	3,4x11,5
					3,8x10	3,4x13
					3,8x11,5	3,4x15

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*The protocol must be customized according to the bone characteristics and the type of surgical procedure to be performed.*

*Very dense bone (D1-D2): it is preferable to use all drills to avoid bone stress.*

*Soft bone (D3-D4): consider under-preparing the site in diameter to achieve greater primary stability.*

Ø 4.25 mm 

MOUNTER COLLO / NECK IL® 1,8 - 2,8 - 3,8 / IL WIDE® 0,9 - 1,8 - 2,8  
TIPO DI BOCCOLA / SLEEVE TYPE 5,50

SEQUENZA PASSAGGI / STEP SEQUENCE	H 6.0 H 8.5 H 10.0 H 11.5 H 13.0 H 15.0					
	1.	2.	3.			
	tissue punch	tissue punch	tissue punch	tissue punch	tissue punch	tissue punch
	crestal mill 4,25	crestal mill 4,25	crestal mill 4,25	crestal mill 4,25	crestal mill 4,25	crestal mill 4,25
	initial drill	initial drill	initial drill	initial drill	initial drill	initial drill
SEQUENZA FRESE / DRILLS SEQUENCE (Ø x H)	2x6	2x6	2x6	2x6	2x6	2x6
	3,4x6	2x8,5	2x8,5	2x8,5	2x8,5	2x8,5
	3,8x6	3,4x8,5	2x10	2x10	2x10	2x10
	4,25x6	3,8x8,5	3,4x10	2x11,5	2x11,5	2x11,5
		4,25x8,5	3,8x10	3,4x10	2x13	2x13
			4,25x10	3,4x11,5	3,4x10	2x15
				3,8x10	3,4x11,5	3,4x10
				3,8x11,5	3,4x13	3,4x11,5
				4,25x10	3,8x10	3,4x13
				4,25x11,5	3,8x11,5	3,4x15
				3,8x13	3,8x10	
				4,25x10	3,8x11,5	
				4,25x11,5	3,8x13	
				4,25x13	3,8x15	
					4,25x10	
					4,25x11,5	
					4,25x13	
					4,25x15	

Ø 5.00 mm 

MOUNTER COLLO / NECK IL® 2,8 - 3,8 / IL WIDE® 0,9 - 1,8 - 2,8  
TIPO DI BOCCOLA / SLEEVE TYPE 5,50

SEQUENZA PASSAGGI / STEP SEQUENCE	H 6.0 H 8.5 H 10.0 H 11.5 H 13.0 H 15.0					
	1.	2.	3.			
	tissue punch	tissue punch	tissue punch	tissue punch	tissue punch	tissue punch
	crestal mill 5,00	crestal mill 5,00	crestal mill 5,00	crestal mill 5,00	crestal mill 5,00	crestal mill 5,00
	initial drill	initial drill	initial drill	initial drill	initial drill	initial drill
SEQUENZA FRESE / DRILLS SEQUENCE (Ø x H)	2x6	2x6	2x6	2x6	2x6	2x6
	3,4x6	2x8,5	2x8,5	2x8,5	2x8,5	2x8,5
	3,8x6	3,4x8,5	2x10	2x10	2x10	2x10
	4,25x6	3,8x8,5	3,4x10	2x11,5	2x11,5	2x11,5
	5x6	4,25x8,5	3,8x10	3,4x10	2x13	2x13
		5x8,5	4,25x10	3,4x11,5	3,4x10	2x15
			5x10	3,8x10	3,4x11,5	3,4x10
				3,8x11,5	3,4x13	3,4x11,5
				4,25x10	3,8x10	3,4x13
				4,25x11,5	3,8x11,5	3,4x15
				5x10	3,8x13	3,8x10
				5x11,5	4,25x10	3,8x11,5
					4,25x11,5	3,8x13
					4,25x13	3,8x15
					5x10	4,25x10
					5x11,5	4,25x11,5
					5x13	4,25x13
						4,25x15
						5x10
						5x11,5
						5x13
						5x15