# ACTIVE

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*Built-in innovation. Trusted clinical performance.*

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**MICROTHREADING**

*Reduces stress on the dental crown and enhances cell adhesion.*

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
**DAS: DOUBLE ACIDIFIED SURFACE**

*Reduced risk of contamination with foreign microparticles. Excellent osseointegration.*

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**TRUNCATED CONE SHAPE**

*Adapts to every clinical need, such as split crest, sinus lift, post-extraction, and immediate loading, thanks to its great primary stability.*

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**VARIABLE GEOMETRY THREAD**

Increases the bone-implant contact surface, increasing the BIC (a Ø 3.8x10 implant develops a surface area of 131 mm<sup>2</sup>, further increased by the 43% roughness index).

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**ROUNDED APEX**

*Prevents damage to important anatomical structures (i.e. Schneiderian membrane or mandibular nerve).*

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*Lanced surgical drill for cortical incision*

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*1.75 mm surgical drill for working length*

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*2.0 mm surgical drill for working length*

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*Ø 3,0 final drill, insert the drill tip completely for 8 mm*

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*Ø 3,8 final drill, insert the drill tip completely for 8 mm*

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*2.8 mm surgical drill for working length*

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*Ø 4,5 final drill, insert the drill tip completely for 8 mm*

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*Ø 5,2 final drill, insert the drill tip completely for 8 mm*

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Ø 6,0 Ø 7,0

*Primarily for post-extraction use, so the protocol will vary based on the clinical situation.*

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
*The Active implant can be used as a tap with the aid of the ratchet, performing three quarters of a turn clockwise and two quarters of a turn counterclockwise.*

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
*Prosthetic components and Instruments*

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

FRESA PILOTA / PILOT DRILL

100152

FRESA PILOTA CON STOP / PILOT DRILL WITH STOP

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*IMPLANT CONNECTION ACTIVE 3.0 HEX FOR RATCHET*

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*IMPLANT CONNECTION ACTIVE 4.0 HEX FOR RATCHET*

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*The indications and images are just general instructions that should be adapted on a case-by-case basis according to the patient's needs. These recommendations are not intended to replace specialized training or the expertise of the clinician, who—also based on personal experience—can make different decisions.*

### **ESSENTIAL PROCEDURES TO ENSURE THE SAFETY, EFFECTIVENESS, AND FUNCTIONALITY OF INSTRUMENTS**

*All surgical instruments must be thoroughly cleaned, disinfected, and sterilized before first use and before each subsequent reuse. It is essential to check the proper functioning of the instruments before each procedure and, if necessary, to immediately replace any instruments showing signs of wear or malfunction. Failure to follow these instructions may increase the risk of infection and intraoperative complications for the patient.*

### **PREPARING THE AREA FOR IMPLANT PLACEMENT**

*Every procedure must be performed gradually and carefully, avoiding sudden movements that could compromise safety. For traditional surgery check the correct positioning of the drill stops, ensuring they are fully tightened. Choose the length in accordance with the implant insertion plan and the given warnings. Periodically check surgical micromotors used with drills and taps to prevent malfunctions. Regularly check the drills for wear and replace them if they lose effectiveness in order to ensure precision and safety during implant site preparation.*

### **PROSTHESIS MAINTENANCE**

*To prevent potential mechanical and biological complications associated with prosthetic implants—which are widely documented in the literature—it is essential for the patient to maintain proper oral hygiene and undergo scheduled periodic checkups. These measures help prolong the functional life of the device. Periodic checks of the calibration of abutment and prosthetic screws are also recommended to ensure the stability and efficiency of the device. If a feeling of instability or malfunction occurs, the patient should be advised to seek prompt medical attention.*

### **DISPOSING OF INSTRUMENTS AND IMPLANTS**

*Surgical instruments, consisting of small metal parts, can be disposed of as metal waste only if properly cleaned. Otherwise, they must be treated as biological waste. Similarly, dental implants removed from the patient's mouth must be treated as biological waste. For each disposal procedure, it is recommended to refer to the local regulations in force regarding the management of medical waste.*