FIVE e TiSmart

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
IMPLANTS WITH UNIQUE CONICAL CONNECTION

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
The implant features a 22° conical connection with an internal hex and is designed for subcrestal placement. According to the surgical protocol, insertion at a depth of 0.5–0.8 mm is recommended. It is distinguished by the design of its threads, which differs between the FIVE and TiSmart lines.

The implant body is cylindrical in the coronal and central portions, while the apical portion has a tapered (truncated-cone) shape. This specific geometry promotes cortical bone healing and assists the surgeon in managing bone spaces, consistently ensuring excellent primary stability and an optimal aesthetic outcome.

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
PLATFORM SWITCHING

The constant connection size, across all diameters, allows for prosthetic switching from 0.3 mm to 1.1 mm, depending on the implant diameter. This design helps reduce peri-implant bone resorption and consequently preserves the soft tissues by respecting the biological width, to the full benefit of long-term aesthetic outcomes.

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
CORONAL MORPHOLOGY

The "Machined" surface in the first 0.5 mm is designed to reduce plaque accumulation, while the constant-pitch microthread significantly contributes to reducing bone resorption.

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
D.S.A. CONNECTION

The Double Seal Action (D.S.A.) connection features a dual-geometry design. The coronal portion is conical, with a 22° angle, acting as a guide for prosthetic insertion and ensuring an optimal bacterial seal between the abutment and the implant. The more internal portion is hexagonal and is intended for abutment stabilization and prosthetic positioning at 60° intervals.

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COILS

FIVE: The Five implant features a single-start, triangular thread with a 0.9 mm pitch. The thread depth increases progressively, ensuring a larger bone-implant contact surface in cancellous bone.

TiSmart: The TiSmart implant has a double-start thread with an asymmetric trapezoidal profile. The variable-geometry threads allow modulation of primary stability in different bone qualities and clinical situations.

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APICAL GROOVES

Three or four grooves, depending on the implant diameter and length, make the implant self-tapping, facilitating insertion. Once osseointegrated, they ensure high anti-rotational stability.

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ROUNDED APEX

The rounded apex helps safeguard anatomical structures at risk, such as the mandibular nerve and the Schneiderian membrane.

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TFRLAN 

11

FRESA CILINDRICA Ø 2,0
CYLINDRICAL DRILL Ø 2,0

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

FIVE & TiSmart Implant Surgical Protocol

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
For illustrative purposes, the surgical sequence for an implant FIVE (TTI) with the following characteristics is shown below:

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
IMPLANT'S SUBCRESTAL EMBEDDING

For correct insertion, it is essential to verify the subcrestal positioning of the implant, with an embedding of 0.8 mm below the bone crest (Fig. 1).

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THE USE OF THE H10 STOP ENSURES A MILLING DEPTH OF 10.8 mm THAT GARANTEES THE INSERTION IMPLANT'S SUBCRESTAL EMBEDDING OF 0,5-0,8 mm. (Fig. 2)

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
TO ACCURATELY DETERMINE THE TOTAL DRILLING DEPTH, IT IS RECOMMENDED TO CONSIDER THE USEFUL HOLE DEPTH PLUS THE LENGTH OF THE DRILL APEX, WHICH VARIES FOR EACH DIAMETER. (Fig. 3)

FIVE - TTI	Ø3,4	Ø3,75	Ø4,1	Ø4,5	Ø5,0
Ti-Smart - TSK	Ø3,5		Ø4,0		Ø4,8

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Length of the drill apex and diameter of the implant associated with each drill

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TSK Prosthetic Components for FIVE & TiSmart

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PER IMPRONTA ANALOGICA FOR ANALOGIC PRINT

TSKTRST

Transfer a strappo
Closed tray transfer

Vite TSKVTRST inclusa / screw TSKVTRST included

TSKTR

Transfer tecnica Pick-up
Pick-up transfer

Vite TSKVTR inclusa / screw TSKVTR included

TSKTRSTSH

Transfer a strappo Short
Short Closed tray transfer

Vite TSKVTRSTSH inclusa / screw TSKVTRSTSH included

TSKTRSH

Transfer tecnica Pick-up Short
Short Pick-up transfer

Vite TSKVTRSH inclusa / screw TSKVTRSH included

MONCONI - ABUTMENTS

MONCONI PROVVISORI TEMPORARY ABUTMENTS

TSKMPN

Moncone provvisorio in titanio non rotante
Temporary titanium abutment, non-rotating

Vite TSKVPR inclusa / screw TSKVPR included


TSKMPR

Moncone provvisorio in titanio rotante
Temporary titanium abutment, rotating

SU RICHIESTA / ON REQUEST: vite TSKVPRTX / screw TSKVPRTX

MONCONI DRITTI STRAIGHT ABUTMENTS

TSKMD15

Moncone dritto h1,5
Straight abutment h1.5

TSKMD30

Moncone dritto h3
Straight abutment h3

Vite TSKVPR inclusa / screw TSKVPR included


TSKMD40

Moncone dritto h4
Straight abutment h4

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
TSKMA25°4


Moncone angolato 25° h4
Angled abutment 25° h4

Vite TSKVPR inclusa / *screw TSKVPR included*

SU RICHIESTA / *ON REQUEST*: vite TSKVPRTX / *screw TSKVPRTX*

T-BASE CROMO COBALTO CON CANNULA CALCINABILE PER PROTESI AVVITATA CHROME COBALT T-BASE WITH CALCINABLE CANNULA FOR SCREW-RETAINED PROSTHESIS

TSKCNCAD

T-Base non rotante in Cromo Cobalto con cannula calcinabile
Cobalt Chrome T-Base with calcinable cannula, non-rotating

TSKCRCAD

T-Base rotante in Cromo Cobalto con cannula calcinabile
Cobalt Chrome T-Base with calcinable cannula, rotating

Vite TSKVPR inclusa / *screw TSKVPR included*

SU RICHIESTA / *ON REQUEST*: vite TSKVPRTX / *screw TSKVPRTX*

COMPONENTI CAD/CAM - CAD/CAM COMPONENTS

SCANBODY INTRAORALE INTRAORAL SCANBODY

TSKIOS

Scanbody
Scanbody

Vite TSKVPR inclusa / screw TSKVPR included

ANALOGO DIGITALE DIGITAL ANALOGUE

TSKANCAD

Analogo digitale
Digital analogue

Vite TISANSCREW inclusa / screw TISANSCREW included

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
The non-rotating bases are all equipped with an anti-rotational lobe. The bases can be cut at the three deepest grooves. Within the library, there are three geometries, each available with four different length offsets.

Vite TSKVPR inclusa / screw TSKVPR included


SU RICHIESTA / ON REQUEST: vite TSKVPRTX / screw TSKVPRTX

BASETTE ROTANTI ROTATING BASES

TSKTIRCAD1

Basetta rotante h1
Rotating base h1

TSKTIRCAD2

Basetta rotante h2
Rotating base h2

Vite TSKVPR inclusa / screw TSKVPR included

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The rotating bases can be cut at the deepest groove. Within the library, there are two geometries, each available with two different length offsets.

Matematiche Digitali disponibili - Available Digital Mathematics

Unico[®]

ABUTMENT

PER LINEA FIVE Ø 3,75/4,1 E LINEA TISMART Ø 4,0
FOR FIVE LINE Ø 3.75/4.1 AND TISMART LINE Ø 4.0

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*Neck 2.8: The Unico Abutment provides a convergent transmucosal section matching the IL[®] implant neck with 2.8 mm height.

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MONCONI DRITTI STRAIGHT ABUTMENTS

UTI3728-6N

Unico Abutment dritto, H6, Collo 2,8*, non rotante
Unico Abutment straight h6, Neck 2.8*, non-rotating

UTI3728-6R

Unico Abutment dritto, H6, Collo 2,8*, rotante
Unico Abutment straight h6, Neck 2.8*, rotating

UTI3728-8N

Unico Abutment dritto, H8, Collo 2,8*, non rotante
Unico Abutment straight h8, Neck 2.8*, non-rotating

UTI3728-8R

Unico Abutment dritto, H8, Collo 2,8*, rotante
Unico Abutment straight h8, Neck 2.8*, rotating

UTI3728-10N

Unico Abutment dritto, H10, Collo 2,8*, non rotante
Unico Abutment straight h10, Neck 2.8*, non-rotating

UTI3728-10R

Unico Abutment dritto, H10, Collo 2,8*, rotante
Unico Abutment straight h10, Neck 2.8*, rotating

MONCONI ANGOLATI ANGLED ABUTMENTS

UTI3715°28-6N

Unico Abutment angolato 15° h6, Collo 2,8*, non rotante
Unico Abutment 15° angled h6, Neck 2.8*, non-rotating

UTI3715°28-6R

Unico Abutment angolato 15° h6, Collo 2,8*, rotante
Unico Abutment 15° angled h6, Neck 2.8*, rotating

PER LINEA FIVE Ø 4,5/5,0 E LINEA TISMART Ø 4,8 FOR FIVE LINE Ø 4.5/5.0 AND TISMART LINE Ø 4.8

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*Neck 2.8: The Unico Abutment provides a convergent transmucosal section matching the IL® implant neck with 2.8 mm height.

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MONCONI DRITTI STRAIGHT ABUTMENTS

UTTI4528-6N

Unico Abutment dritto, H6, Collo 2,8*, non rotante
Unico Abutment straight h6, Neck 2.8*, non-rotating

UTTI4528-6R

Unico Abutment dritto, H6, Collo 2,8*, rotante
Unico Abutment straight h6, Neck 2.8*, rotating

UTTI4528-8N

Unico Abutment dritto, H8, Collo 2,8*, non rotante
Unico Abutment straight h8, Neck 2.8*, non-rotating

UTTI4528-8R

Unico Abutment dritto, H8, Collo 2,8*, rotante
Unico Abutment straight h8, Neck 2.8*, rotating

UTTI4528-10N

Unico Abutment dritto, H10, Collo 2,8*, non rotante
Unico Abutment straight h10, Neck 2.8*, non-rotating

UTTI4528-10R

Unico Abutment dritto, H10, Collo 2,8*, rotante
Unico Abutment straight h10, Neck 2.8*, rotating

MONCONI ANGOLATI ANGLED ABUTMENTS

UTTI4515°28-6N

Unico Abutment angolato 15° h6, Collo 2,8*, non rotante
Unico Abutment 15° angled h6, Neck 2.8*, non-rotating

UTTI4515°28-6R

Unico Abutment angolato 15° h6, Collo 2,8*, rotante
Unico Abutment 15° angled h6, Neck 2.8*, rotating

VITI PROTESICHE - PROSTHETIC SCREW

TSKVTRST

Vite per transfer a strappo
Closed tray transfer screw

TSKVPR

Vite protesica
Prosthetic screw

TSKVTRSTSH

Vite per transfer a strappo Short
Short closed tray transfer screw

TSKVPRTX

Vite protesica torx
Torx prosthetic screw

TSKVTR

Vite per transfer tecnica Pick-up
Pick-up transfer screw

TISANSCREW

Vite per analogo digitale
Screw for digital analogue

TSKVTRSH

Vite per transfer tecnica Pick-up Short
Short Pick-up transfer screw

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The advised prosthetic tightening torque is 25 N-cm

MUA – MONCONI MULTI UNIT - MULTI UNIT ABUTMENT

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
MUA - TRANSFER AND ANALOGUES

TIS900003

Transfer pick up per MUA
Transfer pick up for MUA

TIS900001

Analogo per MUA
Analogue for MUA

Transfer pick up TIS900003: vite TIS900004 inclusa / screw TIS900004 included

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Scanbody TISINTRA: vite TIS100003 inclusa / *screw TIS100003 included*
Analogo CAD / *CAD analogue* TISANCAD: vite TISANSCREW inclusa / *screw TISANSCREW included*

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Moncone provvisorio / *Temporary abutment* TIS100222: vite TIS100003 inclusa / *screw TIS100003 included*

MUA - MONCONE CALCINABILE

MUA - CALCINABLE ABUTMENT

TIS900002

Moncone calcinabile per MUA
Calcinable abutment for MUA

Moncone calcinabile / *Calcinable abutment* TIS900002: vite TIS100003 inclusa / *screw TIS100003 included*

MUA T-BASE CROMO COBALTO CON CANNULA CALCINABILE

COBALT CHROME T-BASE MUA WITH CALCINABLE CANNULA

TISBACCR

T-Base in Cromo Cobalto con cannula
calcinabile per MUA
*Cobalt Chrome T-Base with calcinable cannula
for MUA*

T-Base MUA / *T-Base MUA*: vite TIS100003 inclusa / *screw TIS100003 included*

MUA T-BASE TITANIO

TITANIUM T-BASE MUA

TISBATI

T-Base in titanio per MUA
Titanium T-Base for MUA

TISBATISH

T-Base short in titanio per MUA
Titanium T-Base short for MUA

T-Base MUA / *T-Base MUA*: vite TIS100003 inclusa / *screw TIS100003 included*

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MUA angolato / Angled MUA: vite TSK100001 inclusa / screw TSK100001 included

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Vite per analogo digitale MUA
Screw for MUA digital analogue

TIS900004

Vite transfer per MUA
Transfer screw for MUA

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The advised prosthetic tightening torque is 25 N-cm

TIS100003:

The advised prosthetic tightening torque is 15 N-cm

MONCONI LOCATOR - LOCATOR ABUTMENTS

MONCONI LOCATOR LOCATOR ABUTMENTS

TSKLQ1

Moncone Locator h1
Locator abutment h1

Cappetta ABC/KA-CL-11 inclusa / ABC/KA-CL-11 cap included

TSKLQ2

Moncone Locator h2
Locator abutment h2

Cappetta ABC/KA-CL-11 inclusa / ABC/KA-CL-11 cap included

TSKLQ3

Moncone Locator h3
Locator abutment h3

Cappetta ABC/KA-CL-11 inclusa / ABC/KA-CL-11 cap included

TSKLQ4

Moncone Locator h4
Locator abutment h4

Cappetta ABC/KA-CL-11 inclusa / ABC/KA-CL-11 cap included

TSKLQ5

Moncone Locator h5
Locator abutment h5

Cappetta ABC/KA-CL-11 inclusa / ABC/KA-CL-11 cap included

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TSKMS1

Moncone a sfera h1
Ball abutment h1

Cappetta CN001 inclusa / *CN001 cap included*

TSKMS2

Moncone a sfera h2
Ball abutment h2

Cappetta CN001 inclusa / *CN001 cap included*

TSKMS3

Moncone a sfera h3
Ball abutment h3

Cappetta CN001 inclusa / *CN001 cap included*

TSKMS4

Moncone a sfera h4
Ball abutment h4

Cappetta CN001 inclusa / *CN001 cap included*

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Fresa cilindrica Ø 2,0
Ø 2.0 cylindrical drill

TTISTOP2007 - TTISTOP2085
TTISTOP2010 - TTISTOP2011
TTISTOP2013 - TTISTOP2015

Stop per fresa cilindrica Ø 2,0
Stop for Ø 2.0 cylindrical drill

FRESE FIVE E TISMART FIVE & TISMART DRILLS


TTIFR34

Fresa conica
Ø 3,4 per Five (TTI)
Ø 3,5 per TiSmart (TSK)
*Conical drill
Ø 3.4 for Five (TTI)
Ø 3.5 for TiSmart (TSK)*

TTISTOP3407 - TTISTOP3485
TTISTOP3410 - TTISTOP3411
TTISTOP3413 - TTISTOP3415

Stop per fresa TTIFR34
Stop for TTIFR34 drill

TTIFRC34

Fresa corticale
Ø 3,4 per Five (TTI)
Ø 3,5 per TiSmart (TSK)
*Cortical drill
Ø 3.4 for Five (TTI)
Ø 3.5 for TiSmart (TSK)*

TTIFR37

Fresa conica
Ø 3,75 per Five (TTI)
*Conical drill
Ø 3.75 for Five (TTI)*

TTISTOP3707 - TTISTOP3785
TTISTOP3710 - TTISTOP3711
TTISTOP3713 - TTISTOP3715

Stop per fresa TTIFR37
Stop for TTIFR37 drill

TTIFRC37

Fresa corticale
Ø 3,75 per Five (TTI)
*Cortical drill
Ø 3.75 for Five (TTI)*

TTIFR41

Fresa conica
 Ø 4,1 per Five (TTI)
 Ø 4,0 per TiSmart (TSK)
*Conical drill
 Ø 4,1 for Five (TTI)
 Ø 4,0 for TiSmart (TSK)*

TTISTOP4107 - TTISTOP4185
 TTISTOP4110 - TTISTOP4111
 TTISTOP4113 - TTISTOP4115

Stop per fresa TTIFR41
Stop for TTIFR41 drill

TTIFRC41

Fresa corticale
 Ø 4,1 per Five (TTI)
 Ø 4,0 per TiSmart (TSK)
*Cortical drill
 Ø 4,1 for Five (TTI)
 Ø 4,0 for TiSmart (TSK)*

TTIFR45

Fresa conica
 Ø 4,5 per Five (TTI)
*Conical drill
 Ø 4,5 for Five (TTI)*

TTISTOP4507 - TTISTOP4585
 TTISTOP4510 - TTISTOP4511
 TTISTOP4513 - TTISTOP4515

Stop per fresa TTIFR45
Stop for TTIFR45 drill

TTIFRC45

Fresa corticale
 Ø 4,5 per Five (TTI)
*Cortical drill
 Ø 4,5 for Five (TTI)*


TTIFR50

Fresa conica
 Ø 5,0 per Five (TTI)
 Ø 4,8 per TiSmart (TSK)
*Conical drill
 Ø 5,0 for Five (TTI)
 Ø 4,8 for TiSmart (TSK)*

TTISTOP5007 - TTISTOP5085
 TTISTOP5010 - TTISTOP5011
 TTISTOP5013 - TTISTOP5015

Stop per fresa TTIFR50
Stop for TTIFR50 drill

TTIFRC50

Fresa corticale
 Ø 5,0 per Five (TTI)
 Ø 4,8 per TiSmart (TSK)
*Cortical drill
 Ø 5,0 for Five (TTI)
 Ø 4,8 for TiSmart (TSK)*

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INSTRUMENTS FOR PROSTHETIC COMPONENTS

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TCOTORXC

Cacciavite Torx corto
(non incluso nel kit)
*Short Torx screwdriver
(not included in the kit)*

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The indications and images are just general instructions that should be adapted on a case-by-case basis according to the patient's needs. These recommendations are not intended to replace specialized training or the expertise of the clinician, who—also based on personal experience—can make different decisions.

ESSENTIAL PROCEDURES TO ENSURE THE SAFETY, EFFECTIVENESS, AND FUNCTIONALITY OF INSTRUMENTS

All surgical instruments must be thoroughly cleaned, disinfected, and sterilized before first use and before each subsequent reuse. It is essential to check the proper functioning of the instruments before each procedure and, if necessary, to immediately replace any instruments showing signs of wear or malfunction. Failure to follow these instructions may increase the risk of infection and intraoperative complications for the patient.

PREPARING THE AREA FOR IMPLANT PLACEMENT

Every procedure must be performed gradually and carefully, avoiding sudden movements that could compromise safety. For traditional surgery check the correct positioning of the drill stops, ensuring they are fully tightened. Choose the length in accordance with the implant insertion plan and the given warnings. Periodically check surgical micromotors used with drills and taps to prevent malfunctions. Regularly check the drills for wear and replace them if they lose effectiveness in order to ensure precision and safety during implant site preparation.

PROSTHESIS MAINTENANCE

To prevent potential mechanical and biological complications associated with prosthetic implants—which are widely documented in the literature—it is essential for the patient to maintain proper oral hygiene and undergo scheduled periodic checkups. These measures help prolong the functional life of the device. Periodic checks of the calibration of abutment and prosthetic screws are also recommended to ensure the stability and efficiency of the device. If a feeling of instability or malfunction occurs, the patient should be advised to seek prompt medical attention.

DISPOSING OF INSTRUMENTS AND IMPLANTS

Surgical instruments, consisting of small metal parts, can be disposed of as metal waste only if properly cleaned. Otherwise, they must be treated as biological waste. Similarly, dental implants removed from the patient's mouth must be treated as biological waste. For each disposal procedure, it is recommended to refer to the local regulations in force regarding the management of medical waste.