B.O.P.T.[®]

La Tecnica Verticale

The Vertical Technique

by Ignazio Loi

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The preparation kits for the B.O.P.T.[®] vertical technique by Leader Medica were designed in direct collaboration with Dr. Ignazio Loi and manufactured according to his specific clinical guidelines and operating protocols. Each drill has a specific use in the various stages of preparation up to the creation of the prosthetic abutment.

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The B.O.P.T.[®] is a simplified prosthetic procedure that consists of two fundamental steps:

Space and Form

The technique proposes to eliminate the traditional finishing margin by means of an extra-coronal vertical preparation so that the margin itself is dictated by the future prosthetic restoration, simulating the emergence profile and shape of natural teeth.

The transmucosal compartment that is created after preparation of the sulcus and placement of the provisional restoration defines its finish line and space, which after a period of stabilisation results in a thickening and new tissue conformation: the BIOLOGICAL POLYGON.










Vertical preparation is achieved through the use of dedicated diamond drills of various shapes and grits that define their use in the various preparation steps.

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ASSORTIMENTO FRESE | ASSORTMENT OF DRILLS

	codice code	descrizione description	grana grit	codice colore anello code colour ring
	LM862.105.010	Fresa a fiamma Flame drill	media medium	incolore colourless
	LM862.125.012	Fresa a fiamma Flame drill	grossa coarse	 verde green
	LM862.040.012	Fresa a fiamma Flame drill	fine fine	 rosso red
	LM862.020.012	Fresa a fiamma Flame drill	superfine super fine	 giallo yellow
	LM862.125.016	Fresa a fiamma Flame drill	grossa coarse	 verde green
	LM862.040.016	Fresa a fiamma Flame drill	fine fine	 rosso red
	LM862.020.016	Fresa a fiamma Flame drill	superfine super fine	 giallo yellow
	LM863.125.012	Fresa a fiamma Flame drill	grossa coarse	 verde green
	LM863.040.012	Fresa a fiamma Flame drill	fine fine	 rosso red
	LM863.020.012	Fresa a fiamma Flame drill	superfine super fine	 giallo yellow
	LM863.125.016	Fresa a fiamma Flame drill	grossa coarse	 verde green
	LM863.040.016	Fresa a fiamma Flame drill	fine fine	 rosso red
	LM863.020.016	Fresa a fiamma Flame drill	superfine super fine	 giallo yellow

PREPARATION DRILLS

	LM368.125.023	Fresa a ogiva <i>Ogive drill</i>	grossa <i>coarse</i>	 verde <i>green</i>
	LM368.040.023	Fresa a ogiva <i>Ogive drill</i>	fine <i>fine</i>	 rosso <i>red</i>
	LM368.020.023	Fresa a ogiva <i>Ogive drill</i>	superfine <i>super fine</i>	 giallo <i>yellow</i>
	LM862.105.018	Fresa tronco conica <i>Truncated cone drill</i>	media <i>medium</i>	incolore <i>colourless</i>

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Drills **LM862.125.012**, **LM862.040.012**, **LM862.125.016**, **LM862.040.016**, bear a laser marking to indicate the depth of the preparation.

See the correspondence of each notch on the side.

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Le quote sono espresse in mm. | *The dimensions are quoted in mm.*

Preparation and finishing drills for temporary bridges

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
Drill kit for finishing and polishing temporaries for the B.O.P.T.® technique

The finishing and polishing of the temporary are in the B.O.P.T.® technique the key to clinical and biological success of the technique.

The use of the drills, discs, rubber band and brush in a sequential manner allows their specific use in the various steps up to the final finished and polished temporary restoration. Each drill has a function in each individual step.

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The ball drill is used to unload the temporary within the shell and to remove the various interferences before relining.

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Larger multi-blade drills are used sequentially in the first pass to roughen and remove excess resin from the temporary obtained during relining.

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Multi-blade lanceolate drills are used for finishing the interproximal spaces (large) and for opening papillae spaces (small).

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The spindle-mounted paper disc allows the temporary restoration to be completely finished and pre-polished in all its spaces and surfaces.

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The double-sided diamond disc is specifically for finishing deep interproximal spaces.

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
The pink grit is kneaded with a fine diamond grit that ensures that the final polish achieves a degree of excellence on all surfaces.

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The boar brush used at the end of the various steps, with the use of specific pastes or dry polishing, gives us a final mirror-like polish.

ASSORTIMENTO STRUMENTI | INSTRUMENT ASSORTMENT

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leadermedica.com

Registered office: Via Longhin, 11
35129 Padova Italy

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Leader Medica S.r.l. is responsible for both Certifications and Free Sale Certificates.

The other companies mentioned above are responsible for the world distribution of the devices made by Leader Medica S.r.l.

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LM Medical Technology (Shanghai) Co., Ltd
Sabettimed Srl
Tanto le Certificazioni che i Certificati di Libera Vendita sono in capo a Leader Medica S.r.l.