MITO

IMPIANTO CONICO ESAGONO INTERNO CONICAL IMPLANT WITH INTERNAL HEX CONNECTION

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The MITO implant is made of pure 99% 'cold worked' titanium free of aluminium.

*Four diameters are available:
3.40 - 3.75 - 4.20 - 5.00 with heights ranging from 6 mm to 15 mm.*

Body Implant: The MITO implant has a conical 'core' with dual principle coils with variable geometry; sharper and deeper in the apical part, square-shaped and less pronounced in the more coronal part in contact with the cortical bone. The implant thus maintains a cylindrical profile, which adapts well to all clinical situations, allowing excellent primary stability even in the presence of bone with low density.

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MITO features a 45° angled connection in the most coronal portion, designed to facilitate abutment positioning. The internal hexagon, 2.42 mm wide and 2.2 mm deep, ensures high stability of the prosthetic restoration. The through-screw thread is M 1.8.

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The apical portion features two helical grooves, making the implant self-tapping even in dense bone.

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With a double-lead, self-compacting design, they provide high primary implant stability even in D3/D4 bone. A cutting apical profile transitions into a square geometry in the coronal portion, specifically designed to preserve crestal bone.

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
01

**REDUCED RISK OF
CONTAMINATION WITH
FOREIGN MICROPARTICLES**


02

**EXCELLENT
OSSEOINTEGRATION**

The DAS surface (present on the body of all Leader Medica implants) is a double-acidified surface. This surface is obtained using a chemical method that involves the use of acid mixtures to roughen the titanium surface by subtraction. The first phase involves a surface activation treatment (i.e. the creation of a slight macro-roughness) using hydrofluoric acid. After a thorough washing with demineralized water, the actual etching responsible for the micro-roughness occurs using a mixture of sulfuric and hydrochloric acids. After a further wash, the passivation process takes place: it's a deep cleaning of the surface to restore the passive layer of the material through the controlled supply of oxygen that promotes the formation of titanium oxide. This type of surface allows for an arithmetic mean roughness (Ra) of 1.35 μm . These roughness values are comparable to those expressed by SLA-type surfaces (Sandblasted Largegrid and Acid-etched) and are in perfect agreement with the literature. As demonstrated by Wennerberg A., the long-term results between these two moderately rough surfaces are superimposable. The advantage of the DAS surface (which is much more costly than the SLA one) is that it prevents the deposition of blasting particles on the implant surface, which often results in potentially harmful aluminum oxide residues.

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






We internally certify surface roughness and cleanliness by using scanning electron microscopy (SEM) and X-rays.

The protocol applied to each implant batch was developed with the contribution of Professor Luigi Canullo, recognized as one of the leading experts in the field of implant surfaces and related manufacturing and cleaning processes.

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Industrial processes (turning, surface treatment, packaging) always produce processing residues that deposit on the surface. Although decontamination phases involving washing and ultrasonic baths are included in the various steps of the manufacturing process, microscopic examinations often reveal carbon residues or acid residues. The effect of these contaminants undoubtedly decreases osteoblastic cell adhesion and simultaneously produces an inflammatory response.

The most effective strategy for removing them is to use a biophysical process called plasma. This technology uses an ionized gas or an ionized gas mixture to remove any contaminants from the implant surface.


The plasma used for Leader Medica implants, through specific settings and argon-oxygen gas mixtures, effectively removes organic and inorganic contaminants, making the surface ready for interaction with bone cells.

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Certification of the Protocol for Surface Treatments and Processes

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WHAT IS THE LEADER MEDICA CERTIFIED SURFACE?

The LEADER MEDICA CERTIFIED PROTOCOL concerns the processes and treatments performed on the surfaces of Leader Medica dental implants, as shown in the diagram on the following page.

For the purpose of implant osseointegration, the considered implant surface parameters are, on the one hand, ROUGHNESS, and on the other, CLEANLINESS and PURITY.

Roughness refers to the deviations of the peaks and valleys of the surface asperities calculated along a linear path traced on the surface;

Cleanliness and Purity refer to the absence of localized accumulations of contaminants on the implant surface.

CREATED BY

The LEADER MEDICA CERTIFIED PROTOCOL was developed, refined, tested and validated by the engineers, chemists and technicians of Leader Medica, together with Professor Luigi Canullo (Unicamillus University, Rome).

WHY?

The LEADER MEDICA CERTIFIED PROTOCOL was developed to enhance the quality of surface Roughness, Purity, and Cleanliness, ensuring the predictability and repeatability of these parameters.

IS IT TRUE?

Leader Medica have a latest-generation SEM (Scanning Electron Microscope) equipped with EDS (Energy Dispersive Spectroscopy), used to verify and capture images of the surface Roughness and Cleanliness of dental implants treated according to the LEADER MEDICA CERTIFIED PROTOCOL.

The Roughness images (with 8000x magnification) and Surface Cleanliness images (via energy-dispersive spectrometry) are produced by Leader Medica and always provided to the clinician together with the dental implant belonging to the batch verified by SEM-EDS.


FINAL RESULTS

The results obtained through the processes and treatments of the LEADER MEDICA CERTIFIED PROTOCOL, in terms of Roughness, Purity, and Surface Cleanliness, are considered EXCELLENT:

1. Surface roughness (Ra) measured at 1.3 µm;
2. A characterization of absolute purity and surface cleanliness: (Aluminum: 0.00% and Carbon: 0.04%, based on a verified sample of 5,000 dental implants)

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02 LAVAGGIO CON SOLVENTI
Washing with solvents

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MITO Implant Line

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CONICAL IMPLANT WITH INTERNAL HEX CONNECTION

Ø IMPIANTO	Ø 3.40 mm 	Ø 3.75 mm 	Ø 4.20 mm 	Ø 5.00 mm 
ESAGONO	2,3 mm	*2,42 mm	*2,42 mm	*2,42 mm
Ø APICE	Ø 2,16 mm	Ø 3,15 mm	Ø 3,15 mm	Ø 3,9 mm



H 6.0			MITO4260	MITO5060
H 8.0		MITO3780	MITO4280	MITO5080
H 8.5	MITO3485			
H 10.0	MITO3410	MITO3710	MITO4210	MITO5010
H 11.5	MITO3411	MITO3711	MITO4211	MITO5011
H 13.0	MITO3413	MITO3713	MITO4213	MITO5013
H 15.0	MITO3415	MITO3715	MITO4215	MITO5015

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*Ø 2.42 hex, Screw-Vent-compatible connection

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For illustrative purposes, the surgical sequence for an implant with the following characteristics is shown below:

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
PLACEMENT OF THE SYSTEM AT CRESTAL LEVEL BONE LEVEL

For correct insertion, it is essential to check the crestal Bone Level positioning of the implant, in line with the bone crest (Fig. 1).

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IN ORDER TO DETERMINE THE EXACT MILLING DEPTH, IT IS RECOMMENDED TO CONSIDER THE LENGTH OF THE DRILL TIP, WHICH VARIES FOR EACH DRILL DIAMETER. (Fig. 2)

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TRANSFER E ANALOGHI PER Ø 3,40
TRANSFER AND ANALOGUES FOR Ø 3.40

ILTR
 Transfer pick up
Transfer pick up

Vite ILVTR inclusa / screw ILVTR included

ILANCAD3428
 Analogo uso analogico e CAD
Analogue for conventional and CAD use

Vite ILANSCREW inclusa / screw ILANSCREW included

TRANSFER E ANALOGHI PER Ø 3,75 - 4,20 - 5,00
TRANSFER AND ANALOGUES FOR Ø 3.75 - 4.20 - 5.00

MITOTR
 Transfer pick up
Transfer pick up

Vite MITOVTR inclusa / screw MITOVTR included

MITOTRST
 Transfer a strappo
Closed tray transfer

Vite MITOVPLST inclusa / screw MITOVPLST included

MITOANCAD
 Analogo uso analogico e CAD
Analogue for conventional and CAD use

Vite TISANSCREW128 inclusa / screw TISANSCREW128 included

MONCONI - ABUTMENTS

MONCONI PROVVISORI PER Ø 3,40 TEMPORARY ABUTMENTS FOR Ø 3.40

ILMPN

Moncone provvisorio in titanio non rotante
Temporary titanium abutment, non-rotating

ILMPR

Moncone provvisorio in titanio rotante
Temporary titanium abutment, rotating

Vite ILVPR inclusa / screw ILVPR included

SU RICHIESTA / ON REQUEST: vite ILVPRTX / screw ILVPRTX

MONCONI DRITTI PER Ø 3,40 STRAIGHT ABUTMENTS FOR Ø 3.40

ILMD6

Moncone dritto Profile non rotante h6
Profile straight abutment non-rotating h6

ILMD8

Moncone dritto Profile non rotante h8
Profile straight abutment non-rotating h8

Vite ILVPRSV inclusa / screw ILVPRSV included

MONCONI ANGOLATI PER Ø 3,40 ANGLED ABUTMENTS FOR Ø 3.40

ILMA15°1

Moncone angolato Profile 15° h1
Profile angled abutment 15° h1

ILMA25°1

Moncone angolato Profile 25° h1
Profile angled abutment 25° h1

ILMA15°2

Moncone angolato Profile 15° h2
Profile angled abutment 15° h2

ILMA25°2

Moncone angolato Profile 25° h2
Profile angled abutment 25° h2

Vite ILVPRSV inclusa / screw ILVPRSV included

T-BASE CROMO COBALTO CON CANNULA CALCINABILE PER PROTESI AVVITATA PER Ø 3,40 CHROME COBALT T-BASE WITH CALCINABLE CANNULA FOR SCREW-RETAINED PROSTHESIS FOR Ø 3.40

ILTBCRCON

T-Base non rotante in Cromo Cobalto con cannula calcinabile
Cobalt Chrome T-Base with calcinable cannula, non-rotating

ILTBCRCOR

T-Base rotante in Cromo Cobalto con cannula calcinabile
Cobalt Chrome T-Base with calcinable cannula, rotating

Vite ILVPR inclusa / screw ILVPR included

MONCONI PROVVISORI PER Ø 3,75 - 4,20 - 5,00
TEMPORARY ABUTMENTS FOR Ø 3.75 - 4.20 - 5.00

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The non-rotating bases are all equipped with an anti-rotational lobe. The bases can be cut at the three deepest grooves. Within the library, there are three geometries, each available with four different length offsets.

BASETTE ROTANTI PER Ø 3,75 - 4,20 - 5,00

ROTATING BASES FOR Ø 3.75 - 4.20 - 5.00

MITOTIRCADO

Basetta rotante h0,4
Rotating base h0.4

Vite MITOVPR inclusa / screw MITOVPR included

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Basetta rotante h1,4
Rotating base h1.4

SU RICHIESTA / ON REQUEST: vite MITOVPRTX / screw MITOVPRTX

The rotating bases can be cut at the deepest groove. Within the library, there are two geometries, each available with two different length offsets.

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*Neck 2.8: The Unico Abutment provides a convergent transmucosal section matching the IL[®] implant neck with 2.8 mm height.

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MONCONI DRITTI STRAIGHT ABUTMENTS


UMITO3728-6N

Unico Abutment dritto, H6, Collo 2,8*, non rotante
Unico Abutment straight h6, Neck 2.8*, non-rotating

UMITO3728-6R

Unico Abutment dritto, H6, Collo 2,8*, rotante
Unico Abutment straight h6, Neck 2.8*, rotating

UMITO3728-8N

Unico Abutment dritto, H8, Collo 2,8*, non rotante
Unico Abutment straight h8, Neck 2.8*, non-rotating

UMITO3728-8R

Unico Abutment dritto, H8, Collo 2,8*, rotante
Unico Abutment straight h8, Neck 2.8*, rotating

UMITO3728-10N

Unico Abutment dritto, H10, Collo 2,8*, non rotante
Unico Abutment straight h10, Neck 2.8*, non-rotating

UMITO3728-10R

Unico Abutment dritto, H10, Collo 2,8*, rotante
Unico Abutment straight h10, Neck 2.8*, rotating

MONCONI ANGOLATI ANGLED ABUTMENTS

UMITO3715°28-6N

Unico Abutment angolato 15° h6, Collo 2,8*, non rotante
Unico Abutment 15° angled h6, Neck 2.8*, non-rotating

UMITO3715°28-6R

Unico Abutment angolato 15° h6, Collo 2,8*, rotante
Unico Abutment 15° angled h6, Neck 2.8*, rotating

VITI PROTESICHE - PROSTHETIC SCREW

ILVTR

Vite per transfer pick up Ø 3,4
Pick up transfer Ø 3.4 screw

MITOVTR

Vite per transfer pick up Ø 3,75 - 4,20 - 5,00
Pick up transfer Ø 3.75 - 4.20 - 5.00 screw

MITOVPLST

Vite per transfer a strappo
Closed tray transfer screw

ILANSCREW

Vite per analogo digitale M 1,8
Screw for digital analogue M 1.8

MITOVPR

Vite protesica
Prosthetic screw

MITOVPRTX

Vite protesica torx
Torx prosthetic screw

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The advised prosthetic tightening torque is 25 N-cm

MUA – MONCONI MULTI UNIT - MULTI UNIT ABUTMENT

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
MUA - TRANSFER AND ANALOGUES

IL900003

Transfer pick up per MUA
Transfer pick up for MUA

TIS900001

Analogo per MUA
Analogue for MUA

Transfer pick up IL900003: vite IL900004 inclusa / *screw IL900004 included*

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Scanbody ILINTRA: vite IL100003 inclusa / *screw IL100003 included*

Analogo CAD / *CAD analogue* ILANCADMUA: vite TISANSCREW128 inclusa / *screw TISANSCREW128 included*

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Moncone provvisorio / *Temporary abutment* IL100222: vite IL100003 inclusa / *screw IL100003 included*

MUA - MONCONE CALCINABILE

MUA - CALCINABLE ABUTMENT

IL900002

Moncone calcinabile per MUA
Calcinable abutment for MUA

Moncone calcinabile / *Calcinable abutment* IL900002: vite IL100003 inclusa / *screw IL100003 included*

MUA T-BASE CROMO COBALTO CON CANNULA CALCINABILE

COBALT CHROME T-BASE MUA WITH CALCINABLE CANNULA

ILBACCR

T-Base in Cromo Cobalto con cannula calcinabile per MUA
Cobalt Chrome T-Base with calcinable cannula for MUA

T-Base MUA: vite IL100003 inclusa / *screw IL100003 included*

MUA T-BASE TITANIO

TITANIUM T-BASE MUA

ILBATI

T-Base in titanio per MUA
Titanium T-Base for MUA

ILBATISH

T-Base short in titanio per MUA
Titanium T-Base short for MUA

T-Base MUA: vite IL100003 inclusa / *screw IL100003 included*

VITI PROTESICHE PER MUA - MUA PROSTHETIC SCREW

IL100003

Vite protesica per MUA
Prosthetic screw for MUA

TIS100003TX

Vite protesica torx per MUA
Torx prosthetic screw for MUA

MITOVPR

Vite protesica per MUA angolato
Prosthetic screw for angled MUA

IL900004

Vite transfer per MUA
Transfer screw for MUA

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The advised prosthetic tightening torque is 25 N-cm

IL100003:

The advised prosthetic tightening torque is 15 N-cm

MONCONI LOCATOR - LOCATOR ABUTMENTS

MONCONI LOCATOR PER Ø 3,40

LOCATOR ABUTMENTS FOR Ø 3.40

ILLO1

Moncone Locator h1
Locator abutment h1

ILLO2

Moncone Locator h2
Locator abutment h2

ILLO3

Moncone Locator h3
Locator abutment h3

ILLO4

Moncone Locator h4
Locator abutment h4

ILLO5

Moncone Locator h5
Locator abutment h5

MONCONI LOCATOR PER Ø 3,75 - 4,20 - 5,00

LOCATOR ABUTMENTS FOR Ø 3.75 - 4.20 - 5.00

MITOLO1

Moncone Locator h1
Locator abutment h1

Cappetta ABC/KA-CL-11 inclusa / ABC/KA-CL-11 cap included

MITOLO2

Moncone Locator h2
Locator abutment h2

Cappetta ABC/KA-CL-11 inclusa / ABC/KA-CL-11 cap included

MITOLO3

Moncone Locator h3
Locator abutment h3

Cappetta ABC/KA-CL-11 inclusa / ABC/KA-CL-11 cap included

MITOLO4

Moncone Locator h4
Locator abutment h4

Cappetta ABC/KA-CL-11 inclusa / ABC/KA-CL-11 cap included

MITOLO5

Moncone Locator h5
Locator abutment h5

Cappetta ABC/KA-CL-11 inclusa / ABC/KA-CL-11 cap included

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MITOMS1

Moncone a sfera h1
Ball abutment h1

Cappetta CN001 inclusa / CN001 cap included

CA011-50

Porta cappetta in metallo
Metal cap holder

MITOMS2

Moncone a sfera h2
Ball abutment h2

Cappetta CN001 inclusa / CN001 cap included

CN001

Cappetta rosa Ø 2,5 per moncone a sfera
Pink cap Ø 2,5 for ball abutment

MITOMS3

Moncone a sfera h3
Ball abutment h3

Cappetta CN001 inclusa / CN001 cap included

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Fresa cilindrica Ø 2,0
Ø 2.0 cylindrical drill

TSTOPC2006 - TSTOPC2080
TSTOPC2010 - TSTOPC2011
TSTOPC2013 - TSTOPC2015

Stop per fresa cilindrica Ø 2,0
Stop for Ø 2.0 cylindrical drill

FRESE MITO
MITO DRILLS

ILFR34

Fresa conica Ø 3,40
Ø 3.40 Conical drill

TSTOPC2806 - TSKSTOP3585
TSKSTOP3510 - TSKSTOP3511
TSKSTOP3513 - TSKSTOP3515

Stop per fresa conica Ø 3,40
Stop for Ø 3.40 conical drill

ILFRC34

Fresa corticale Ø 3,40
Ø 3.40 Cortical drill

MITOFR37

Fresa conica Ø 3,75
Ø 3.75 Conical drill

MITOSTOP3706 - MITOSTOP3780
MITOSTOP3710 - MITOSTOP3711
MITOSTOP3713 - MITOSTOP3715

Stop per fresa conica Ø 3,75
Stop for Ø 3.75 conical drill

MITOFRC37

Fresa corticale Ø 3,75
Ø 3.75 Cortical drill

MITOFR42

Fresa conica Ø 4,20
Ø 4.20 Conical drill

MITOSTOP4206 - MITOSTOP4280
MITOSTOP4210 - MITOSTOP4211
MITOSTOP4213 - MITOSTOP4215

Stop per fresa conica Ø 4,20
Stop for Ø 4.20 conical drill

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INSTRUMENTS FOR PROSTHETIC COMPONENTS

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The indications and images are just general instructions that should be adapted on a case-by-case basis according to the patient's needs. These recommendations are not intended to replace specialized training or the expertise of the clinician, who—also based on personal experience—can make different decisions.

ESSENTIAL PROCEDURES TO ENSURE THE SAFETY, EFFECTIVENESS, AND FUNCTIONALITY OF INSTRUMENTS

All surgical instruments must be thoroughly cleaned, disinfected, and sterilized before first use and before each subsequent reuse. It is essential to check the proper functioning of the instruments before each procedure and, if necessary, to immediately replace any instruments showing signs of wear or malfunction. Failure to follow these instructions may increase the risk of infection and intraoperative complications for the patient.

PREPARING THE AREA FOR IMPLANT PLACEMENT

Every procedure must be performed gradually and carefully, avoiding sudden movements that could compromise safety. For traditional surgery check the correct positioning of the drill stops, ensuring they are fully tightened. Choose the length in accordance with the implant insertion plan and the given warnings. Periodically check surgical micromotors used with drills and taps to prevent malfunctions. Regularly check the drills for wear and replace them if they lose effectiveness in order to ensure precision and safety during implant site preparation.

PROSTHESIS MAINTENANCE

To prevent potential mechanical and biological complications associated with prosthetic implants—which are widely documented in the literature—it is essential for the patient to maintain proper oral hygiene and undergo scheduled periodic checkups. These measures help prolong the functional life of the device. Periodic checks of the calibration of abutment and prosthetic screws are also recommended to ensure the stability and efficiency of the device. If a feeling of instability or malfunction occurs, the patient should be advised to seek prompt medical attention.

DISPOSING OF INSTRUMENTS AND IMPLANTS

Surgical instruments, consisting of small metal parts, can be disposed of as metal waste only if properly cleaned. Otherwise, they must be treated as biological waste. Similarly, dental implants removed from the patient's mouth must be treated as biological waste. For each disposal procedure, it is recommended to refer to the local regulations in force regarding the management of medical waste.

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Intellectual Property

B.O.P.T.[®]

HPA[®]

BODYBLADE[®]

FMC[®]

CONO-IN[®]

3P
SMART
DEVICES

LEADER MEDICA[®]
曆 德

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Certifications and Authorisations

EU-MDD
Certified

EU-MDR
Certified

Ente certificante | Certification Authority

UNI CEI EN ISO 13485:2021

Ente certificante | Certification Authority