
MITO

MADE IN ITALY

Comitato Scientifico & Opinion Leader

/ SCIENTIFIC COMMITTEE & KEY OPINION LEADER

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The DNA of Evolution

Leader Medica pursues a Project based on the following core concepts:

- **Maximizing Quality** in every phase of the dental implant devices' designing and manufacturing;
- **Training** for anyone who interacts with the company (i.e. internal collaborators, sales associates, dentists, dental technicians, etc.);
- **Digitalization** of both the processes involved in the production of medical devices and the processes relevant to the use of the devices (e.g., guided surgery systems);
- Generating and/or strengthening a close **Relationship of Trust** between the Company and the Client with the aim of creating the foundation for a successful and synergistic business and professional development, characterized by valuable training and digital content.

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Device quality, training, digitalization, and internationalization are the key concepts of the evolution that sets us apart and defines us for the future.

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Leader Medica has been designing and manufacturing high-precision implant systems and prosthetic components since 1994, overseeing every phase in-house, from the research to the production as well as the clinical and training support. Innovation is part of our DNA: we study, discover, and create solutions that simplify and improve therapeutic efficacy. We are a dynamic, constantly evolving company focused on international trade and boasting a network of outstanding scientific collaborations.

LEADER MEDICA consolidates its growth

In 2024, Leader Medica acquired 3P Implafavourite, a company that has distinguished itself in the market with the 3P Smart Devices brand and as the creator of the CONO-IN® implant line featuring a conical Morse-taper connection, a true expression of technological and design excellence. Through this acquisition, a second manufacturing facility is established in Pinerolo (Turin, Italy), strengthening the Group's industrial capacity and strategic footprint.

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
01 REDUCED RISK OF
CONTAMINATION WITH
FOREIGN MICROPARTICLES

02 EXCELLENT
OSSEOINTEGRATION

The DAS surface (present on the body of all Leader Medica implants) is a double-acidified surface. This surface is obtained using a chemical method that involves the use of acid mixtures to roughen the titanium surface by subtraction. The first phase involves a surface activation treatment (i.e. the creation of a slight macro-roughness) using hydrofluoric acid. After a thorough washing with demineralized water, the actual etching responsible for the micro-roughness occurs using a mixture of sulfuric and hydrochloric acids. After a further wash, the passivation process takes place: it's a deep cleaning of the surface to restore the passive layer of the material through the controlled supply of oxygen that promotes the formation of titanium oxide. This type of surface allows for an arithmetic mean roughness (Ra) of 1.35 μm . These roughness values are comparable to those expressed by SLA-type surfaces (Sandblasted Largegrid and Acid-etched) and are in perfect agreement with the literature. As demonstrated by Wennerberg A., the long-term results between these two moderately rough surfaces are superimposable. The advantage of the DAS surface (which is much more costly than the SLA one) is that it prevents the deposition of blasting particles on the implant surface, which often results in potentially harmful aluminum oxide residues.

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






MICRORUGOSITÀ - MICROROUGHNESS

SEM | Mag= 8.000 X Signal A= SE1

We internally certify surface roughness and cleanliness by using scanning electron microscopy (SEM) and X-rays.

The protocol applied to each implant batch was developed with the contribution of Professor Luigi Canullo, recognized as one of the leading experts in the field of implant surfaces and related manufacturing and cleaning processes.

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
Industrial processes (turning, surface treatment, packaging) always produce processing residues that deposit on the surface. Although decontamination phases involving washing and ultrasonic baths are included in the various steps of the manufacturing process, microscopic examinations often reveal carbon residues or acid residues. The effect of these contaminants undoubtedly decreases osteoblastic cell adhesion and simultaneously produces an inflammatory response.

The most effective strategy for removing them is to use a biophysical process called plasma. This technology uses an ionized gas or an ionized gas mixture to remove any contaminants from the implant surface.


The plasma used for Leader Medica implants, through specific settings and argon-oxygen gas mixtures, effectively removes organic and inorganic contaminants, making the surface ready for interaction with bone cells.

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Certification of the Protocol for Surface Treatments and Processes

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WHAT IS THE LEADER MEDICA CERTIFIED SURFACE?

The LEADER MEDICA CERTIFIED PROTOCOL concerns the processes and treatments performed on the surfaces of Leader Medica dental implants, as shown in the diagram on the following page.

For the purpose of implant osseointegration, the considered implant surface parameters are, on the one hand, ROUGHNESS, and on the other, CLEANLINESS and PURITY.

Roughness refers to the deviations of the peaks and valleys of the surface asperities calculated along a linear path traced on the surface;

Cleanliness and Purity refer to the absence of localized accumulations of contaminants on the implant surface.

CREATED BY

The LEADER MEDICA CERTIFIED PROTOCOL was developed, refined, tested and validated by the engineers, chemists and technicians of Leader Medica, together with Professor Luigi Canullo (Unicamillus University, Rome).

WHY?

The LEADER MEDICA CERTIFIED PROTOCOL was developed to enhance the quality of surface Roughness, Purity, and Cleanliness, ensuring the predictability and repeatability of these parameters.

IS IT TRUE?

Leader Medica have a latest-generation SEM (Scanning Electron Microscope) equipped with EDS (Energy Dispersive Spectroscopy), used to verify and capture images of the surface Roughness and Cleanliness of dental implants treated according to the LEADER MEDICA CERTIFIED PROTOCOL.

The Roughness images (with 8000x magnification) and Surface Cleanliness images (via energy-dispersive spectrometry) are produced by Leader Medica and always provided to the clinician together with the dental implant belonging to the batch verified by SEM-EDS.


FINAL RESULTS

The results obtained through the processes and treatments of the LEADER MEDICA CERTIFIED PROTOCOL, in terms of Roughness, Purity, and Surface Cleanliness, are considered EXCELLENT:

1. Surface roughness (Ra) measured at 1.3 µm;
2. A characterization of absolute purity and surface cleanliness: (Aluminum: 0.00% and Carbon: 0.90%, based on a verified sample of 5,000 dental implants)

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02 LAVAGGIO CON SOLVENTI
Washing with solvents

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
The MITO implant features a conical "core" with variable-geometry, double-start threads; sharper and deeper at the apical portion, square-shaped and less pronounced at the coronal portion in contact with the cortical bone.

This design maintains a cylindrical profile, making it suitable for all clinical situations and ensuring excellent primary stability even in low-density bone.

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CONNECTION

MITO features a 90° angled connection in the most coronal portion, designed to facilitate abutment positioning. The internal hexagon, 2.42 mm wide and 2.2 mm deep, ensures high stability of the prosthetic restoration. The through-screw thread is M 1.8.

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THREADS

With a double-lead, self-compacting design, they provide high primary implant stability even in D3/D4 bone. A cutting apical profile transitions into a square geometry in the coronal portion, specifically designed to preserve crestal bone.

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
APICAL NOTCHES

The apical portion features two helical grooves, making the implant self-tapping even in dense bone.

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DAS SURFACE - DOUBLE ACIDIFIED SURFACE

The implant root surface undergoes a double acid-etching treatment designed to optimize wettability and adhesiveness and to promote optimal osseointegration.

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MITO Implant Line

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CONICAL IMPLANT WITH INTERNAL HEX CONNECTION

Ø IMPIANTO Ø IMPLANT	Ø 3.40 mm 	Ø 3.75 mm 	Ø 4.20 mm 	Ø 5.00 mm 
ESAGONO HEXAGON	2,3 mm	*2,42 mm	*2,42 mm	*2,42 mm
Ø APICE Ø APEX	Ø 2,16 mm	Ø 3,15 mm	Ø 3,15 mm	Ø 3,9 mm


H 6.0			MITO4260	MITO5060
H 8.0		MITO3780	MITO4280	MITO5080
H 8.5	MITO3485			
H 10.0	MITO3410	MITO3710	MITO4210	MITO5010
H 11.5	MITO3411	MITO3711	MITO4211	MITO5011
H 13.0	MITO3413	MITO3713	MITO4213	MITO5013
H 15.0	MITO3415	MITO3715	MITO4215	MITO5015

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**Ø 2.42 hex, Screw-Vent-compatible connection*

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For illustrative purposes, the surgical sequence for an implant with the following characteristics is shown below:

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
PLACEMENT OF THE SYSTEM AT CRESTAL *BONE LEVEL*

For correct insertion, it is essential to check the crestal Bone Level positioning of the implant, in line with the bone crest (Fig. 1).

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IN ORDER TO DETERMINE THE EXACT MILLING DEPTH, IT IS RECOMMENDED TO CONSIDER THE LENGTH OF THE DRILL TIP, WHICH VARIES FOR EACH DRILL DIAMETER. (Fig. 2)

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*Neck 2.8: The Unico Abutment provides a convergent transmucosal section matching the IL[®] implant neck with 2.8 mm height.

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MONCONI DRITTI STRAIGHT ABUTMENTS

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ILTR

Transfer pick up
Transfer pick up

Vite ILVTR inclusa / screw ILVTR included

ILANCAD3428

Analogo uso analogico e CAD
Analogue for analogic and CAD use

Vite ILANScrew inclusa / screw ILANScrew included

TRANSFER E ANALOGHI PER Ø 3,75 - 4,20 - 5,00

TRANSFER AND ANALOGUES FOR Ø 3.75 - 4.20 - 5.00

MITOTR

Transfer pick up
Transfer pick up

Vite MITOVTR inclusa / screw MITOVTR included

MITOTRST

Transfer a strappo
Closed tray transfer

Vite MITOVPLST inclusa / screw MITOVPLST included

MITOANCAD

Analogo uso analogico e CAD
Analogue for analogic and CAD use

Vite TISANScrew128 inclusa / screw TISANScrew128 included

MONCONI - ABUTMENTS

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TEMPORARY ABUTMENTS FOR Ø 3.40

ILMPN

Moncone provvisorio in titanio non rotante
Temporary titanium abutment, non-rotating

ILMPR

Moncone provvisorio in titanio rotante
Temporary titanium abutment, rotating

Vite ILVPR inclusa / screw ILVPR included

SU RICHIESTA / ON REQUEST: vite ILVPRTX / screw ILVPRTX

MONCONI DRITTI PER Ø 3,40

STRAIGHT ABUTMENTS FOR Ø 3.40

ILMD6

Moncone dritto Profile non rotante h6
Profile straight abutment non-rotating h6

ILMD8

Moncone dritto Profile non rotante h8
Profile straight abutment non-rotating h8

Vite ILVPRSV inclusa / screw ILVPRSV included

MONCONI ANGOLATI PER Ø 3,40

ANGLED ABUTMENTS FOR Ø 3.40

ILMA15°1

Moncone angolato Profile 15° h1
Profile 15° angled abutment h1

ILMA25°1

Moncone angolato Profile 25° h1
Profile 25° angled abutment h1

ILMA15°2

Moncone angolato Profile 15° h2
Profile 15° angled abutment h2

ILMA25°2

Moncone angolato Profile 25° h2
Profile 25° angled abutment h2

Vite ILVPRSV inclusa / screw ILVPRSV included

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COBALT-CHROME T-BASE WITH CALCINABLE CANNULA FOR SCREW-RETAINED PROSTHESIS FOR Ø 3.40

ILTBCRCON

T-Base non rotante in Cromo Cobalto con cannula calcinabile
Cobalt-Chrome T-Base with calcinable cannula, non-rotating

ILTBCRCOR

T-Base rotante in Cromo Cobalto con cannula calcinabile
Cobalt-Chrome T-Base with calcinable cannula, rotating

Vite ILVPR inclusa / screw ILVPR included

MONCONI PROVVISORI PER Ø 3,75 - 4,20 - 5,00

TEMPORARY ABUTMENTS FOR Ø 3.75 - 4.20 - 5.00

MITOMPN

Moncone provvisorio in titanio non rotante
Temporary titanium abutment, non-rotating

MITOMPR

Moncone provvisorio in titanio rotante
Temporary titanium abutment, rotating

Vite MITOVPR inclusa / screw MITOVPR included

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COBALT-CHROME T-BASE WITH CALCINABLE CANNULA FOR SCREW-RETAINED PROSTHESIS FOR Ø 3.75 - 4.20 - 5.00

MITOCROCO

T-Base non rotante in Cromo Cobalto con cannula calcinabile
Cobalt-Chrome T-Base with calcinable cannula, non-rotating

MITOCROCOR

T-Base rotante in Cromo Cobalto con cannula calcinabile
Cobalt-Chrome T-Base with calcinable cannula, rotating

Vite MITOVPR inclusa / screw MITOVPR included

COMPONENTI CAD/CAM - CAD/CAM COMPONENTS

SCANBODY PER Ø 3,40

SCANBODY FOR Ø 3.40

ILIOS

Scanbody

Scanbody

Vite ILVPR inclusa / screw ILVPR included

ANALOGO DIGITALE PER Ø 3,40

DIGITAL ANALOGUE FOR Ø 3.40

ILANCAD3428

Analogo uso analogico e CAD

Analogue for analogic and CAD use

Vite ILANSCREW inclusa / screw ILANSCREW included

BASETTE NON ROTANTI PER Ø 3,40

NON-ROTATING BASES FOR Ø 3.40

ILTB1N-4

Basetta non rotante h 1,5

Non-rotating base h1.5

Vite ILVPR inclusa / screw ILVPR included

ILTB2N-4

Basetta non rotante h 2,5

Non-rotating base h2.5

SU RICHIESTA / ON REQUEST: vite ILVPRTX / screw ILVPRTX

BASETTE ROTANTI PER Ø 3,40

ROTATING BASES FOR Ø 3.40

ILTB1R-4

Basetta rotante h 1,5

Rotating base h1.5

Vite ILVPR inclusa / screw ILVPR included

ILTB2R-4

Basetta rotante h 2,5

Rotating base h2.5

SU RICHIESTA / ON REQUEST: vite ILVPRTX / screw ILVPRTX

SCANBODY PER Ø 3,75 - 4,20 - 5,00

SCANBODY FOR Ø 3.75 - 4.20 - 5.00

MITOIOS

Scanbody
Scanbody

Vite MITOVPR inclusa / screw MITOVPR included

ANALOGO DIGITALE PER Ø 3,75 - 4,20 - 5,00

DIGITAL ANALOGUE FOR Ø 3.75 - 4.20 - 5.00

MITOANCAD

Analogo uso analogico e CAD
Analogue for analogic and CAD use

Vite TISANSCREW128 inclusa / screw TISANSCREW128 included

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The non-rotating bases are all equipped with an anti-rotational lobe. The bases can be cut at the three deepest grooves. Within the library, there are three geometries, each available with four different length offsets.

BASETTE ROTANTI PER Ø 3,75 - 4,20 - 5,00

ROTATING BASES FOR Ø 3.75 - 4.20 - 5.00

MITOTIRCADO

Basetta rotante h0,4
Rotating base h0.4

Vite MITOVPR inclusa / *screw MITOVPR included*

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Basetta rotante h1,4
Rotating base h1.4

SU RICHIESTA / *ON REQUEST*: vite MITOVPRTX / *screw MITOVPRTX*

The rotating bases can be cut at the deepest groove. Within the library, there are two geometries, each available with two different length offsets.

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ILVPRTX

Vite protesica torx testa cilindrica per sigillo piano
Torx cylindrical head prosthetic screw for flat seal

TISANSCREW128

Vite per analogo digitale
Screw for digital analogue

ILVTR

Vite per transfer pick up Ø 3,4
Pick up transfer Ø 3.4 screw

ILANSCREW

Vite per analogo digitale M 1,8
Screw for digital analogue M 1.8

MITOVTR

Vite per transfer pick up Ø 3,75 - 4,20 - 5,00
Pick up transfer Ø 3.75 - 4.20 - 5.00 screw

MITOVPR

Vite protesica
Prosthetic screw

MITOVPLST

Vite per transfer a strappo
Closed tray transfer screw

MITOVPRTX

Vite protesica torx
Torx prosthetic screw

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The advised prosthetic tightening torque is 25 N-cm

IL100003:

The advised prosthetic tightening torque is 15 N-cm

MONCONI LOCATOR - LOCATOR ABUTMENTS

MONCONI LOCATOR PER Ø 3,40

LOCATOR ABUTMENTS FOR Ø 3.40

ILLO1

Moncone Locator h1
Locator abutment h1

ILLO2

Moncone Locator h2
Locator abutment h2

ILLO3

Moncone Locator h3
Locator abutment h3

ILLO4

Moncone Locator h4
Locator abutment h4

ILLO5

Moncone Locator h5
Locator abutment h5

MONCONI LOCATOR PER Ø 3,75 - 4,20 - 5,00

LOCATOR ABUTMENTS FOR Ø 3.75 - 4.20 - 5.00

MITOLO1

Moncone Locator h1
Locator abutment h1

Cappetta ABC/KA-CL-11 inclusa / ABC/KA-CL-11 cap included

MITOLO2

Moncone Locator h2
Locator abutment h2

Cappetta ABC/KA-CL-11 inclusa / ABC/KA-CL-11 cap included

MITOLO3

Moncone Locator h3
Locator abutment h3

Cappetta ABC/KA-CL-11 inclusa / ABC/KA-CL-11 cap included

MITOLO4

Moncone Locator h4
Locator abutment h4

Cappetta ABC/KA-CL-11 inclusa / ABC/KA-CL-11 cap included

MITOLO5

Moncone Locator h5
Locator abutment h5

Cappetta ABC/KA-CL-11 inclusa / ABC/KA-CL-11 cap included

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MITOMS1

Moncone a sfera h1
Ball abutment h1

Cappetta CN001 inclusa / CN001 cap included

MITOMS2

Moncone a sfera h2
Ball abutment h2

Cappetta CN001 inclusa / CN001 cap included

MITOMS3

Moncone a sfera h3
Ball abutment h3

Cappetta CN001 inclusa / CN001 cap included

MITOMS4

Moncone a sfera h4
Ball abutment h4

Cappetta CN001 inclusa / CN001 cap included

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Fresa cilindrica Ø 2,0
Ø 2.0 cylindrical drill

TSTOPC2006 - TSTOPC2080
TSTOPC2010 - TSTOPC2011
TSTOPC2013 - TSTOPC2015

Stop per fresa cilindrica Ø 2,0
Stop for Ø 2.0 cylindrical drill

FRESE MITO

MITO DRILLS

ILFR34

Fresa conica Ø 3,40
Ø 3.40 Conical drill

TSTOPC2806 - TSKSTOP3585
TSKSTOP3510 - TSKSTOP3511
TSKSTOP3513 - TSKSTOP3515

Stop per fresa conica Ø 3,40
Stop for Ø 3.40 conical drill

ILFRC34

Fresa corticale Ø 3,40
Ø 3.40 Cortical drill

MITOFR37

Fresa conica Ø 3,75
Ø 3.75 Conical drill

MITOSTOP3706 - MITOSTOP3780
MITOSTOP3710 - MITOSTOP3711
MITOSTOP3713 - MITOSTOP3715

Stop per fresa conica Ø 3,75
Stop for Ø 3.75 conical drill

MITOFRC37

Fresa corticale Ø 3,75
Ø 3.75 Cortical drill

MITOFR42

Fresa conica Ø 4,20
Ø 4.20 Conical drill

MITOSTOP4206 - MITOSTOP4280
MITOSTOP4210 - MITOSTOP4211
MITOSTOP4213 - MITOSTOP4215

Stop per fresa conica Ø 4,20
Stop for Ø 4.20 conical drill

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INSTRUMENTS FOR PROSTHETIC COMPONENTS

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Juxta-sinusal zygomatic implantology, an evolution of the intra-sinusal technique with palatal access, allows simplification of the surgical procedure, reduction of intraoperative time, limitation of postoperative complications, and improvement of prosthetic function. In specific and well-defined clinical situations, the zygomatic implant becomes the only therapeutic solution for the patient.

Leader Medica has therefore developed the MITO Zygomatic system, which incorporates this innovation and includes a zygomatic implant with differentiated surfaces designed to simplify surgery, as well as a surgical kit specifically developed for the juxta-sinusal technique, containing diamond drills used to create the bone "slot" that guides the final drill. A complete prosthetic component range, including MUA abutments up to 60°, allows prosthetic management of all different implant emergence profiles.

The implant features a completely cylindrical morphology and is available in two diameters: 3.75 mm, suitable in cases of reduced zygomatic bone volume, and 4.20 mm, to be used when sufficient zygomatic bone is present. Leader Medica provides dental practitioners with a comprehensive theoretical and practical training program aimed at deepening knowledge of this surgical procedure and enabling its selection as a therapeutic indication in immediate-load implant-prosthetic rehabilitation of patients with severe maxillary atrophy.

Ø esagono 2,42 compatibile connessione Screw-Vent
Ø 2.42 hex, Screw-Vent-compatible connection

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A cylindrical zone of variable height with a smooth, machined surface. It is minimally invasive, with soft tissues in direct contact. Two diameters are available: Ø 3.75 mm and Ø 4.20 mm.

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
Threaded zone for anchorage in the zygomatic bone, 18 mm in height. Double acid-etched surface treatment, with a double-lead thread and a self-tapping profile. It allows easy and stable placement in the zygomatic bone.

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

The MITO Zygomatic surgical kit, available for 3.75 mm and 4.20 mm diameters, is made of autoclavable and therefore sterilizable materials. The innovative drills included in the MITO Zygomatic surgical kit are the result of an extensive study carried out over the last four years by a team of clinical, scientific, and mechanical engineering consultants. This applies to both the steel drills and the diamond drills, whose surface treatments have been specifically developed for the MITO Zygomatic implant system. The drills allow implant site preparation through extremely precise and rapid consecutive steps, ensuring high effectiveness and, above all, strict compliance with protocol guidelines, thereby guaranteeing that the bone does not undergo any alteration or damage. The surgical kit includes all the surgical instruments required for the placement of zygomatic implants and prosthetic components. The manual screwdriver has been studied, designed, and manufactured by Leader Medica specifically for the MITO Zygomatic implant system, ensuring easy and high-performance use of the entire zygomatic system, even in the most complex surgical situations.

ZYTRAY375

Kit chirurgico impianto zigomatico
Zygomatic Implant Surgical Kit
Ø 3,75

ZYTRAY420

Kit chirurgico impianto zigomatico
Zygomatic Implant Surgical Kit
Ø 4,20

MITO Zigomatico, Nasale e Pterigoideo

Componenti Protetiche

MITO Zygomatic, Nasal and Pterygoid, Prosthetic Components

MUA – MONCONI MULTI UNIT - MULTI UNIT ABUTMENT

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
Moncone provvisorio / *Temporary abutment* IL100222: vite IL100003 inclusa / *screw IL100003 included*

MUA - MONCONE CALCINABILE

MUA - CALCINABLE ABUTMENT

IL900002

Moncone calcinabile per MUA
Calcinable abutment for MUA

Moncone calcinabile / *Calcinable abutment* IL900002: vite IL100003 inclusa / *screw IL100003 included*

Matematiche Digitali disponibili - *Available Digital Mathematics*

MUA T-BASE IN CROMO COBALTO CON CANNULA CALCINABILE

COBALT-CHROME T-BASE MUA WITH CALCINABLE CANNULA

ILBACCR

T-Base in Cromo Cobalto con cannula calcinabile per MUA
Cobalt-Chrome T-Base with calcinable cannula for MUA

T-Base MUA: vite IL100003 inclusa / screw IL100003 included

MUA T-BASE TITANIO

TITANIUM T-BASE MUA

ILBATI

T-Base in titanio per MUA
Titanium T-Base for MUA

T-Base MUA: vite IL100003 inclusa / screw IL100003 included

ILBATISH

T-Base short in titanio per MUA
Titanium T-Base short for MUA

SU RICHIESTA / ON REQUEST:
vite TIS100003TX / screw TIS100003TX

VITI PROTESICHE PER MUA - MUA PROSTHETIC SCREW

IL100003

Vite protesica per MUA
Prosthetic screw for MUA

TIS100003TX

Vite protesica torx per MUA
Torx prosthetic screw for MUA

MITOVPR

Vite protesica per MUA angolato
Prosthetic screw for angled MUA

IL900004

Vite transfer per MUA
Transfer screw for MUA

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The indications and images are just general instructions that should be adapted on a case-by-case basis according to the patient's needs. These recommendations are not intended to replace specialized training or the expertise of the clinician, who—also based on personal experience—can make different decisions.

ESSENTIAL PROCEDURES TO ENSURE THE SAFETY, EFFECTIVENESS, AND FUNCTIONALITY OF INSTRUMENTS

All surgical instruments must be thoroughly cleaned, disinfected, and sterilized before first use and before each subsequent reuse. It is essential to check the proper functioning of the instruments before each procedure and, if necessary, to immediately replace any instruments showing signs of wear or malfunction. Failure to follow these instructions may increase the risk of infection and intraoperative complications for the patient.

PREPARING THE AREA FOR IMPLANT PLACEMENT

Every procedure must be performed gradually and carefully, avoiding sudden movements that could compromise safety. For traditional surgery check the correct positioning of the drill stops, ensuring they are fully tightened. Choose the length in accordance with the implant insertion plan and the given warnings. Periodically check surgical micromotors used with drills and taps to prevent malfunctions. Regularly check the drills for wear and replace them if they lose effectiveness in order to ensure precision and safety during implant site preparation.

PROSTHESIS MAINTENANCE

To prevent potential mechanical and biological complications associated with prosthetic implants—which are widely documented in the literature—it is essential for the patient to maintain proper oral hygiene and undergo scheduled periodic checkups. These measures help prolong the functional life of the device. Periodic checks of the calibration of abutment and prosthetic screws are also recommended to ensure the stability and efficiency of the device. If a feeling of instability or malfunction occurs, the patient should be advised to seek prompt medical attention.

DISPOSING OF INSTRUMENTS AND IMPLANTS

Surgical instruments, consisting of small metal parts, can be disposed of as metal waste only if properly cleaned. Otherwise, they must be treated as biological waste. Similarly, dental implants removed from the patient's mouth must be treated as biological waste. For each disposal procedure, it is recommended to refer to the local regulations in force regarding the management of medical waste.

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Certifications and Authorisations

EU-MDD
Certified

EU-MDR
Certified

Ente certificante | Certification Authority

Ente certificante | Certification Authority

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leadermedica.com

Registered office: Via Longhin, 11
35129 Padova Italy

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Leader Medica S.r.l. is responsible for both Certifications and Free Sale Certificates.

The other companies mentioned above are responsible for the world distribution of the devices made by Leader Medica S.r.l.

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LM Medical Technology (Shanghai) Co., Ltd
Sabettimed Srl
Tanto le Certificazioni che i Certificati di Libera Vendita sono in capo a Leader Medica S.r.l.