


# Linea implantare CONO-IN®


## CONO-IN® Implant Line

GUIDA ALLA LETTURA DEGLI IMPIANTI  
GUIDE TO READING IMPLANTS

- H** ALTEZZA IMPIANTO  
IMPLANT HEIGHT
- ØA** DIAMETRO IMPIANTO  
IMPLANT DIAMETER
- ØB** DIAMETRO APICALE  
APICAL DIAMETER



Ø A 3,0 mm ● 3,4 mm ● 3,8 mm ● 4,5 mm ● 5,0 mm ● 6,2 mm ●



Ø B 1,5 mm 1,5 mm 1,6 mm 1,8 mm 2,55 mm 2,65 mm

H 6,5 mm				C45006	C50006	C62006
H 8 mm	C30008	C34008	C38008	C45008	C50008	C62008
H 10 mm	C30010	C34010	C38010	C45010	C50010	C62010
H 12 mm	C30012	C34012	C38012	C45012	C50012	C62012
H 14 mm	C30014	C34014	C38014	C45014	C50014	
H 16 mm	C30016	C34016	C38016	C45016		
H 18 mm		C34018	C38018			
H 20 mm		C34020	C38020			

### VITI TAPPO H 1,5 CAP SCREWS H 1.5

Ø 3.0



Ø3,0  
500258

Ø≥3.4



SMALL  
500261

Ø≥3.4



MEDIUM  
500262

Ø 3.0

Per/for Cono-In® Ø 3.0

Ø≥3.4

Per/for Cono-In®  
Ø 3.4 - Ø 3.8 - Ø 4.5  
Ø 5.0 - Ø 6.2