









Linea Implantare

ACTIVE

ACTIVE Implant Line

- H** ALTEZZA IMPIANTO
IMPLANT HEIGHT
- ØA** DIAMETRO IMPIANTO
IMPLANT DIAMETER
- ØB** DIAMETRO APICALE
APICAL DIAMETER



	Ø A	2,5 mm	3,3 mm	3,8 mm	4,5 mm	5,2 mm	6,0 mm	7,0 mm
								
	Ø B	0,60 mm	1,50 mm	1,60 mm	1,80 mm	2,55 mm	2,80 mm	3,65 mm
H 6 mm				A38006	A45006	A52006	A60006	
H 8 mm			A33008	A38008	A45008	A52008	A60008	A70008
H 10 mm		A25010	A33010	A38010	A45010	A52010	A60010	
H 12 mm		A25012	A33012	A38012	A45012	A52012	A60012	
H 14 mm		A25014	A33014	A38014	A45014			
H 16 mm			A33016	A38016				

MONCONI DA SALDATURA WELDING ABUTMENTS



ES-HEX 3
100095







ES-HEX 4
100096

Impianti ACTIVE ESTETICO

ACTIVE ESTETICO Implants



- H** ALTEZZA IMPIANTO
IMPLANT HEIGHT
- ØA** DIAMETRO IMPIANTO
IMPLANT DIAMETER
- ØB** DIAMETRO APICALE
APICAL DIAMETER




	Ø A	3,8 mm	4,5 mm	5,2 mm
				
	Ø B	1,60 mm	1,80 mm	2,55 mm
H 8 mm				AE52008
H 10 mm		AE38010	AE45010	AE52010
H 12 mm		AE38012	AE45012	
H 14 mm		AE38014	AE45014	

Impianti ORING MICRO


ORING Implants

	$\varnothing A$ 2,5 mm	3,3 mm
		
	$\varnothing B$ 0,65 mm	1,50 mm
	$\varnothing C$ 1,80 mm	1,80 mm
H 10 mm	025010	033010
H 12 mm	025012	033012
H 14 mm	025014	033014



Impianti ORING NORMO

ORING NORMO Implants

	$\varnothing A$ 3,3 mm
	
	$\varnothing B$ 1,50 mm
	$\varnothing C$ 2,50 mm
H 10 mm	ON33010
H 12 mm	ON33012
H 14 mm	ON33014

