

# Linea implantare FIVE

## *FIVE Implant Line*

**IMPIANTI CONNESSIONE CONICA UNICA**  
*IMPLANTS WITH UNIQUE CONICAL CONNECTION*

<b>Ø IMPIANTO</b>	Ø 3.40 mm 	Ø 3.75 mm 	Ø 4,10 mm 	Ø 4,50 mm 	Ø 5.00 mm 
<b>Ø APICE</b>	Ø 2.7 mm	Ø 3.1 mm	Ø 3.4 mm	Ø 3.8 mm	Ø 4.3 mm



<b>H 7.0</b>		TTI3707L	TTI4107L	TTI4507L	TTI5007L
<b>H 8.5</b>	TTI3485L	TTI3785L	TTI4185L	TTI4585L	TTI5085L
<b>H 10.0</b>	TTI3410L	TTI3710L	TTI4110L	TTI4510L	TTI5010L
<b>H 11.5</b>	TTI3411L	TTI3711L	TTI4111L	TTI4511L	TTI5011L
<b>H 13.0</b>	TTI3413L	TTI3713L	TTI4113L	TTI4513L	TTI5013L
<b>H 15.0</b>	TTI3415L	TTI3715L	TTI4115L	TTI4515L	TTI5015L

# Linea implantare TiSmart

## *TiSmart Implant Line*

**IMPIANTI CONNESSIONE CONICA UNICA**  
*IMPLANTS WITH UNIQUE CONICAL CONNECTION*

Ø IMPIANTO	Ø 3.50 mm 	Ø 4,00 mm 	Ø 4.80 mm 
Ø APICE	Ø 2.8 mm	Ø 3.0 mm	Ø 3.8 mm



H 7,0	TSK3507L	TSK4007L	TSK4807L
H 8,5	TSK3585L	TSK4085L	TSK4885L
H 10	TSK3510L	TSK4010L	TSK4810L
H 11,5	TSK3511L	TSK4011L	TSK4811L
H 13,0	TSK3513L	TSK4013L	TSK4813L
H 15,0	TSK3515L	TSK4015L	TSK4815L