



Impianto IL[®] Zigomatico, Nasale e Pterigoideo

IL[®] Zygomatic, Nasal and Pterygoid Implant



IL[®] ZIGOMATICO Ø 3,75 IL[®] ZYGOMATIC Ø 3.75

h 35 mm	ILZY3735	ILZY4235
h 37,5 mm	ILZY3737	ILZY4237
h 40 mm	ILZY3740	ILZY4240
h 42,5 mm	ILZY3742	ILZY4242
h 45 mm	ILZY3745	ILZY4245
h 47,5 mm	ILZY3747	ILZY4247
h 50 mm	ILZY3750	ILZY4250
h 52,5 mm	ILZY3752	ILZY4252
h 55 mm	ILZY3755	ILZY4255
h 60 mm	ILZY3760	ILZY4260



IL[®] NASALE Ø 4,00 IL[®] NASAL Ø 4.00

h 24 mm	ILNS4024
h 26 mm	ILNS4026
h 28 mm	ILNS4028
h 30 mm	ILNS4030



IL[®] PTERIGOIDEO Ø 4,00 IL[®] PTERYGOID Ø 4.00

h 18 mm	ILPT4018
h 20 mm	ILPT4020
h 22 mm	ILPT4022

IL[®] Zigomatico, Nasale e Pterigoideo

Componenti Protetiche

IL[®] Zygomatic, Nasal and Pterygoid, Prosthetic Components

MUA – MONCONI MULTI UNIT - MULTI UNIT ABUTMENT

MUA ANGOLATO ANGLED MUA



ILWMUA17°3

MUA angolato 17° h3 IL WIDE®
17° angled MUA h3 IL WIDE®



ILWMUA17°5

MUA angolato 17° h5 IL WIDE®
17° angled MUA h5 IL WIDE®



ILWMUA30°3

MUA angolato 30° h3 IL WIDE®
30° angled MUA h3 IL WIDE®



ILWMUA30°5

MUA angolato 30° h5 IL WIDE®
30° angled MUA h5 IL WIDE®



ILWMUA45°5

MUA angolato 45° h5 IL WIDE®
Angled MUA 45° h5 IL WIDE®



ILWMUA52°5

MUA angolato 52° h5 IL WIDE®
Angled MUA 52° h5 IL WIDE®



ILWMUA60°5

MUA angolato 60° h5 IL WIDE®
Angled MUA 60° h5 IL WIDE®



TPOMFA

Inseritore per MUA angolato
Angled MUA inserter

MUA angolato / Angled MUA: vite ILVPRMUA inclusa / screw ILVPRMUA included



MUA - VITI DI GUARIGIONE

MUA - HEALING SCREWS



IL100304

Vite di guarigione per MUA h4
Healing screw for MUA h4



IL100308

Vite di guarigione per MUA h8
Healing screw for MUA h8



IL100306

Vite di guarigione per MUA h6
Healing screw for MUA h6

MUA - TRANSFER E ANALOGHI

MUA - TRANSFER AND ANALOGUES



IL900003

Transfer pick up per MUA
Transfer pick up for MUA



TIS900001

Analogo per MUA
Analogue for MUA

Transfer pick up IL900003: vite IL900004 inclusa / *screw IL900004 included*

MUA - SCANBODY E ANALOGHI

MUA - SCANBODY AND ANALOGUES



ILINTRA

Scanbody per MUA
Scanbody for MUA



ILANCADMUA

Analogo CAD per MUA
CAD analogue for MUA

Scanbody ILINTRA: vite IL100003 inclusa / *screw IL100003 included*
Analogo CAD / *CAD analogue* ILANCADMUA: vite TISANSCREW128 inclusa / *screw TISANSCREW128 included*

MUA - MONCONE PROVVISORIO

MUA - TEMPORARY ABUTMENT



IL100222

Moncone provvisorio per MUA
Temporary abutment for MUA

Moncone provvisorio / *Temporary abutment* IL100222: vite IL100003 inclusa / *screw IL100003 included*

MUA - MONCONE CALCINABILE

MUA - CALCINABLE ABUTMENT



IL900002

Moncone calcinabile per MUA
Calcinable abutment for MUA

Moncone calcinabile / *Calcinable abutment* IL900002: vite IL100003 inclusa / *screw IL100003 included*



Matematiche Digitali disponibili - *Available Digital Mathematics*

MUA T-BASE IN CROMO COBALTO CON CANNULA CALCINABILE

COBALT-CHROME T-BASE MUA WITH CALCINABLE CANNULA



ILBACCR

T-Base in Cromo Cobalto con cannula calcinabile per MUA
Cobalt-Chrome T-Base with calcinable cannula for MUA

T-Base MUA: vite IL100003 inclusa / screw IL100003 included

MUA T-BASE TITANIO

TITANIUM T-BASE MUA



ILBATI

T-Base in titanio per MUA
Titanium T-Base for MUA

T-Base MUA: vite IL100003 inclusa / screw IL100003 included



ILBATISH

T-Base short in titanio per MUA
Titanium T-Base short for MUA



SU RICHIESTA / ON REQUEST:
vite TIS100003TX / screw TIS100003TX

VITI PROTESICHE PER MUA - MUA PROSTHETIC SCREW



IL100003

Vite protesica per MUA
Prosthetic screw for MUA



TIS100003TX

Vite protesica torx per MUA
Torx prosthetic screw for MUA



ILVPRMUA

Vite esagonale per MUA angolato
Hexagonal screw for angled MUA



IL900004

Vite transfer per MUA
Transfer screw for MUA



TISANSCREW128

Vite per analogo digitale MUA
Screw for MUA digital analogue

Si consiglia coppia di serraggio protesica a 25 N-cm

IL100003:

Si consiglia coppia di serraggio protesica a 15 N-cm

The advised prosthetic tightening torque is 25 N-cm

IL100003:

The advised prosthetic tightening torque is 15 N-cm